

Product datasheet for **MR231297**

Ifi203 (NM_001302649) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ifi203 (NM_001302649) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ifi203
Synonyms:	Ifi203a; Ifi203b; Mndal
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide
Sequence:**

>MR231297 representing NM_001302649
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGATGGCTGAATACAAGAATATTGTTCTTCTAAAAGGATTAGAGAACATGGAGGACTATCAGTTTAGGA
 CAGTCAAGTCCTTACTAAGAAAAGAACTAAAACCTTACTAAAAAATGCAAGAAGATTATGATAGGATTCA
 GCTGGCGGACTGGATTGAGGACAAGTCCCCAAAGATGCTGGACTAGACAAACTGATAAAAGTGTGTGAG
 CACATAAAAGATCTTAAAGACCTTGCTAAAAAAGTAAAAACAGAGAAGGCAAAAGTTCAAAAAGAAAAGG
 AAGGAAAATGCAAAACGGCTGGGAAAAAAGGGCAAGATGAACTCAGTTCTTCTGAATCTTTTTCAA
 CAATAAAGAAATCAACAAGAGTGTACCGTCTTCAAAGAAAAGAGAAAACAGATCACAAAACTGAAGGT
 GGAAGAAAAGAAAGCTCACCCAGGAGCAGGCTCAGCTCCAGAACTTCAGGAACCAACATAAAAAAG
 AGGAAGATTGCCTCCAGAATCCTCATAAGTCTCCACCAACACCATCCAGCAGTTCTCAAACAAGAAGAA
 TCCTGTTGTGGCGTCTTCTTGGTGGCAAGTCGGGCGTCTACAAGTGGTCCAAAACCCGGGAACAAGAT
 ATCTTAAGGCACAAACAAGACCCTCAAGAAAGGGAGGATTCAAAACTACATTTTGTGAGAGAATTAAGA
 AACAGAAAAACAGAACCATCCGAAAACATAGTATCTCAAAGACAGAGGGTCTCAGGAAAAGCAACAACG
 TGTAGAATTGTCATCAATCAGCAACTTCCAAGCTGTAAGCGAACTCCTGACCTTTGAAGGACTTTTCAGCA
 ATTCTTCGAGCAGACTTCAGTCTTCTCAGAAGCCTCTAGAAGCTCACTTGGATTTAAAGATGTCTCCAA
 TCTCTCCAACATCACCTTGCCACAATCTCCAGTGTCTCTATCCTCAGACTCAAATGTTTCATCTGAACTC
 TAATGCACATTCATCCTGTTGAGTGGTGGCCAGGTTCTCATGTACCATCAGCAACAGGATTTCAGTAAT
 GTCGGGGTCCCTCTTATGCCTTCAGAAAAGTGTTCAGTCTTCCAGTGTCTCCTCAGATATCTCAAATAA
 CAGTACCTAGCAGCATCCGGGACCTTTACCTTCTAAACCAGCAGCATTCAATAGCACAGGCTCCTCA
 CAAACAGGCAATAGCATCCAGAAAATGGCCAACTACTAGAGCGCCATCAGCAAAACTGACAAGCAATCC
 CAGTTCCTAAAGTTGCATCCTCCAGCAACAGCATCCAGCAATTCTCAGGTTCTCATATTCTAGCAACCA
 TGCCAGAAAACATCTCTGTACACAGGTACCTCAAACAAGAACATCCAGTAGCATTTCAGACTCTTAACTC
 TGCTACAGTAAAAGCATCCAAGAATGTCCAGGCCCTCAAGTATCTTTACCAATACGACCAATTATTTT
 CTGGCTTCTGTGGCACCTCCATCAGCAACACCCAAACCGTCTTACATCTCATGGGCCACTTTTGTCAATAG
 CCACAAGCAGAGCCAGACTACTCAAATGAACCAGGAACAGAATCCATTTGTGTCCAGGCTCTTCATGC
 TCCTCCATCAACAGTGTCCAGAAAATAAGTGCACCACTCCATTGACACAAGGGGCAGAATCCAGTACTGGG
 GAAGCCCTTCCCCTTCTAAAGTAAAAGCATGCACGATATCCAAGCACCTCGGGTGTCTTCACCAACTG
 CATCCAGCAGTATCATGAAGCCACATGCAACTTCATCAAAAACATCCAGCAGTCTCCTGGATCAACATGC
 AACTTCATCAACAGCATCCAGCAGCCGTCTTCTGGCTCTAAAGTTGTCTCCAGTAACAGCATCTAGGGCT
 CTCAGTGTATTCCAGGGCTTTTCAGAAAACATACACAGAAGTCCATCCAGCGCACCAAGGAGAGGAACTG
 TACCAAAGGAGCCTTCTAGAGAAGAAGGTCAACATCAAGGTCCCAACAAGTGTGGTGTGAAAGTAAAC
 AGAACCATTTACATATGACTTTGAAGAAACGAAAAGGATGTTTCATGCTACAGTGGCTACTGAAACTGAG
 TTCTTCAGAGTGAAGGTTTTTGACACAGCTCTAATGAGCAAGTTCATCCAGGAAAAGATCATTGCCATAT
 CACATTAATTTGGGTGCAATGGGTTTCTGGAGATATACAGAGCTTCTGTGTGTCGATGTGAACATTA
 CCCAACAATGATTATCTCAAATACACTGAGTGAAGTGTCTATTGCAACTCCAAAAATTTCTATCTTCTC
 TCACAAGCAAAGGGGACATTTGTAATGGAGAATTTGTAGTATTTAAGAAAAGTGAAGGATGAGTGCA
 TTTGCTATGGAATTGGAGATGATACAGGAAAAATGGCAGTGGTGGTTTATGGACGACTCACCAATGTCAG
 GTGTGAACCAGGCAGTAACTTAGACTTGTCTGCTTTGAATTGACTTCCACTAAAGATGTGTGCCTCCTG
 AGGTCTGTAAGGCACAGTTACATGCAGGTCAATGAAGGAAAGCCACTCAACCAGACTCAGTCAGGA
 GAAACTCTCTGGAACCATACTTC

ACGCGTACGCGGCCGCTCGAGCAGAAAACATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR231297 representing NM_001302649
 Red=Cloning site Green=Tags(s)


MMAEYKNIVLLKGLNEMEDYQFRTVKSLLRKELKLTKKMQEDYDRIQLADWIEDKFKPDAGLDKLIKVCE
 HIKDLKDLAKKLTKEKAVQKKKEGKCKTAGKKKGQDELSSESLFNNKESNKSVPSSKKKRQITKTEG
 GKKKKL TQEQAQLPETS G T N I K K E E D C L Q N P H K S P P T P S S S S N K N P V V A S F L L A S R A S T S G A K T R E Q D
 I L R H K Q R P S R K G G F K T T F C E R I K K Q K N R T I R K H S I S K T E G P Q E K Q Q R V E L S S I S N F Q A V S E L L T F E G L S A
 I P S S R L Q S S Q K P L E A H L D L K M S P I S P T S P C H N L P V S L S D S N V H L N S N A H S I L F S G A Q V P H V P S A T G F S N
 V G V P L M P S E T V F S S S S A P Q I S Q I T V P S S I R D L Y L P K P A A F N S T E A P H K Q A I A S R N G Q T T R A P S A K L T S K S
 Q F L K L H P P A T A S S N S Q V P H I L A T M P R N I S V T Q V P Q T R T S S S I Q T L N S A T V K A S K N V Q A P Q V S L P I R P N Y F
 L A S V A P P S A T P N R L T S H G P L L S I A T S R A Q T T Q I E P G T E S I C V Q A L H A P P S T V S R N K C T T P L T Q G A E S S T G
 E A L P L P K V K A C T I F Q A P R V S S P T A S S I M K P H A T S S K T S S S L L D Q H A T S S T A S S R L L A K L S P V T A S R A
 L S A I P G L S E T I H R S P S S A P R R G T V P K E P S R E E G H H Q G P K Q V M V L K V T E P F T Y D F E E T K R M F H A T V A T E T E
 F F R V K V F D T A L M S K F I P G K I I A I S H Y I G C N G F L E I Y R A S C V S D V N I N P T M I I S N T L S E S A I A T P K I S Y L L
 S Q A K G T F V N G E F V V F K K S E R H E C I C Y G I G D D T G K M A V V V Y G R L T N V R C E P G S K L R L V C F E L T S T K D V C L L
 R S V R H S Y M Q V I N E G K P L N P D S V R R N S L E P Y F

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001302649

ORF Size: 2613 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001302649.1](#), [NP_001289578.1](#)

RefSeq Size: 4883 bp

RefSeq ORF: 2616 bp

Locus ID: 15950

UniProt ID: [O35368](#)

Cytogenetics: 1 80.76 cM

MW: 95.8 kDa