

## Product datasheet for **MR231322**

### Uspl1 (NM\_001286825) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Uspl1 (NM_001286825) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Uspl1
Synonyms:	E430026A01Rik
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>MR231322 representing NM\_001286825  
 Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGCCAGCTAAGAGCAGATGTTTACCACCTTTGCCAGACTCTGTGTGCCAGTGGAAAAACCCAGGCC  
 TGTGCTGGTTAGACTGTATCCTGTCAGCCCTGGTGCACCTGGAGGTGTTGAGGAAGACTGTGCTGGAGGC  
 GTGCTCCAGGGAGGAGTGTGATTTGGAAGGCTGTTTGAATGTACCACCAAGCAGATGAGCTGCTGCAC  
 ACCCATCACCTGCATGGCGTCACAGGTGAAGATTGTAATAAATTGACATCAGAAATATTTACAGAAATAG  
 ATACCTGTTTGAATAAAGTTAGAGATGAAATTTTGTAACTTCAACCGAAGCTTAGATGCACATTAGG  
 TGACATGGAAAGTCTGTGTTTGCACCTCCTGTACTGTTAAAGCTTGAACCCATGTTGAAAGCCTCTTT  
 ACATATCTTTTTCTTGGAAATTTGAATGTTCCATTGTGGACACCAGTACCAAAACAGGTGTGTGAAGA  
 GTCTGGTCACCTTACCAATATTGTTCCGTGAGTGGCATCCACTCAATGCTGCCATTTTGGTCCATGTAA  
 CAGCTGCAACAGTAAATCACAAATAAGAAAAATGGTGTGGAAAGAGCGTCGCCATACTCATGCTGCAC  
 TTCGTGGAGGGCTTACCAAGGAGGGATTGCAACACTATGCCTTCCACTTTGAAGGAAGCCTTTATCAAG  
 TCACGTCCGTAATACAGTACCAAGCAAATAACCATTTTATAACATGGTTTTTAGATGCTGATGGGAGTTG  
 GCTGGAGTGTGATGACCTAAAGGGTCCCTGTGCTAAAAGGCACGTGACCTGTGAAGTGCCAGCTTCAGAG  
 ACACATATCGTTATTTGGGAAAGAAAAATCCAGGTGCCAATTGAAGAAGCTGCCTGCCTCCATGCATGA  
 AGCCAAATGTGCAGCCTGTATCAGGTGAGGAGCAGCCGACTTGTCCAGCACTGTGTTCTCTGGCTGGGAC  
 CGCCACCTCAGAACCTCAGTCGCCATCCACATCTATGGCAGGCGCTCCTCAGACTCTCCAGAAATC  
 CAAGCTGTAGCTCATGGAGATAGTGTACTTTCAGGAGCAAAGGCATGGTCGACAGCATTTTACCCTCAG  
 CACTTGAAGAAACCATCCAGGAACTGCCTCTGTTTCTCAGGTTGACTCCAAAGATTGCCTGTTGGAAG  
 CAAACCCGTGGCAGGAAGTGCAGCACTGGTCAGAGTACTTGTCTTTTTCAGCCACAGGACTCCCCGGGTCT  
 TCGGGGTCTTCTTAGTGTCCAGCCTGTGTGAAGGAAAGCTAGTTGCCCATGTGTGGATTCAAGTTTCC  
 CATCCCAAGCTGTGAGCACAGACCTGCAGGCAGTACTGTCCAGGCAGGGGACACTGTGGTTCTTAACCC  
 TGTGACTGATGCCCTGTTCTGTTCTTGTTCAGGAATTGAAGTCCCTGGCAACTGAGAAGGACTCTCAG  
 ACACAGTTGCTGCCACTGAAAACGAGAAGTTAGATCCTGAGCAACCTGGTAAATCTCAGGCATCTAATT  
 TGAGAAAACGAGAAACACAGCCTTCTAAAACAGTAGCAGCCAGGTGACGACAGAATCAGCCTCGGAA  
 AGAAGATCAGAAGAGAGCCTTCGTGGGAAGTTGGGTGAAGGGGCTATTGAGCAGGGGTGGTCTTTTATG  
 CCAACCTGTGTTTGTCTCAGAGTAGAGCAGTAAGTGACCTGCAACCTTCTGTTAAAGGAGCCAGTAATT  
 TTGATGGCTTTAAAACGAAAAGTATAAGCCGAGGTCTAAAAGAATGTCCAGGAAAGCTAAGCACATGGA  
 GAACTGTCTCCAAGGAACAGCTCTCCACCATTGAGCTGGACAGCTGCTCTGACCCAGGCTGCTGAGAAC  
 GCCACCTCAGCTCTGTTGAGGGAGCAGGAAGGCTCACGTCCAGCTCCCCTCAGGCACAGATCTCCTGGCA  
 ATGAAAGCGCCATCTCTCCAGCAAGCCGTGGAGATGCGGCTGAAGATCAGGTTCTAAGCTTCGGCTCAA  
 ACTTCTTAAAAAATTAAGGCAAAAAAGAAAGAACTAGCTGCCCTCATTCTTCCCCGCACCGTGAACCA  
 TCTCTAAGTGACCACTCAGAGCCTGCGTCCCATTGTGGAACACCTGCAAGTGACCACTCAGAGCCTGTAT  
 CCCATTGTGGATCTCAAATGACTGTGAGTCCATAGAAGATCTGTTAAAAGAGCTACAGCATCAAATTGA  
 CTTTGCAGACAGCAAATCTGGGTGCACCTCTGCTCCTGATGCCACCTCAAATAACAGCCAGAGTCATGAG  
 GAGATTTTAGCAGAGTTGCTGTGCCTACAGCCATGTGAGAGCCCTCAGAAAGTGGGAGCTTGAGTTGA  
 GGTACTTGGAAATGGGGGACAGCACCCAGCACACACCTAGTGAATTCAGTGTGCTCTCCAGAACAC  
 ATGTCTGAAACAGGACCATGATTACTGCAGCCCTGAGAAAGGCCAGCGGGAAGTTGATCTGCATTGAGTC  
 ATGGACAGTGCCTGTATAAGAACCTTGAACCTGGGGAGTCCCATGAAGACTGATATCTTTGATGATTTTT  
 TTTCCACCTCAGCACTGAATCTTTAACAAACGACACATTAGACATACCTCATTTTGTGACTCTCTGTT  
 TGAGAATTGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



## Plasmid Map:



ACCN: NM\_001286825

ORF Size: 2670 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001286825.1](#), [NP\\_001273754.1](#)

RefSeq Size: 3714 bp

RefSeq ORF: 2673 bp

Locus ID: 231915

UniProt ID: [Q3ULM6](#)

Cytogenetics: 5 G3

**MW:** 97.7 kDa

**Gene Summary:** SUMO-specific isopeptidase involved in protein desumoylation. Specifically binds SUMO proteins with a higher affinity for SUMO2 and SUMO3 which it cleaves more efficiently. Also able to process full-length SUMO proteins to their mature forms (By similarity). Plays a key role in RNA polymerase-II-mediated snRNA transcription in the Cajal bodies (By similarity). Is a component of complexes that can bind to U snRNA genes (By similarity).[UniProtKB/Swiss-Prot Function]