

## Product datasheet for **MR231344**

### Snap91 (NM\_001277986) Mouse Tagged ORF Clone

#### Product data:

|                    |  |
|--------------------|--|
| Product Type:      | Expression Plasmids                          |
| Product Name:      | Snap91 (NM_001277986) Mouse Tagged ORF Clone |
| Tag:               | Myc-DDK                                      |
| Symbol:            | Snap91                                       |
| Synonyms:          | 91kDa; AP180; F1-20; mKIAA0656               |
| Vector:            | pCMV6-Entry (PS100001)                       |
| E. coli Selection: | Kanamycin (25 ug/mL)                         |
| Cell Selection:    | Neomycin                                     |



[View online »](#)

ORF Nucleotide  
Sequence:


>MR231344 representing NM\_001277986  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGTCGGGCCAAACGCTCACGGATCGGATCGCCGCCGCTCAGTACAGCGTACTGGCTCTGCTGTAGCAA  
GAGCAGTCTGCAAAGCCACCACTCATGAGGTGATGGGCCCAAGAAGAAGCACCTGGACTATTTGATCCA  
GGTACCAATGAGACCAATGTCAATATCCCTCAGATGGCCGACACCCTCTTTGAGCGGGCAGCAAACAGT  
AGCTGGGTGGTGGTATTTAAAGCTTTAGTGACCACACACCATCTCATGGTGCATGGAAATGAGAGATTTA  
TTCAGTATTTGGCCTCTAGGAATACGCTATTCAATCTCAGCAACTTTCTAGATAAAAAGTGGATCCCACGG  
TTATGATATGTCTACGTTACACGGCGTTACAGTAGATACTTGAATGAAAAGGCTTTCTCCTACAGACAG  
ATGGCATTGACTTTGCCAGAGTGAAGAAAGGGGCTGACGGTGCATGAGGACGATGGTTCCTGAAAAGC  
TCCTGAAGAGTATGCCAATCCTGCAGGGGCAGATCGATGCACTGCTGGAGTTTGTATGTCATCCAAATGA  
ACTAACCAATGGTGTCAAATGCTGCATTTATGCTTCTTTTCAAAGATCTTATCAAAGTGGTTCCTGCTGC  
TACAATGACGGCGTCATTAACCTTACTTGAAAAATTTTTCGAGATGAAGAAAGGTCAATGCAAAGACGCGC  
TAGAAATTTACAAGCGATTTCTAACTAGAATGACGAGGGTGTCTGAATTCCTCAAGGTCGCCGAGCAAGT  
TGGTATTGATAAAGGTGACATTTCCGACCTCACACAGGCTCCCAGCAGTCTTATGGAGACCCTTGAACAA  
CATCTAAATACCTTAGAAGGAAAGAAACCTGGAACAATGAAGGATCTGGTGTCCCTCTCCACTAAGTA  
AGTCTTCTCCAGCCAACTGTTACATCTCTAATTCTACACCAGCTAAAATATCGACACATCCCCGCC  
AGTTGACATATTTGCAACAGCATCCGCGGCTGCCCCAGTCAGCTCTGCTAAGCCATCAAGCGATCTCCT  
GATCTTCAGCCGACTTCTCTGGAGCAGCTGCGGGGCAGCAGCACCTGTAGTGCCTCCTCTGGGGGTG  
CGACCCCTGGGAGACCTTTTGGGAGAGGATTCCTTGGTGCATTTCCCTGTTCCTGTGGAAGCACC  
GATTTTCAGACCCATTTGCACCAGAGCCTTCCCCTCCTACTACAACCACTGACCTGCTTCAGCCTCTGCC  
TCGACCACACAGCTGTGACGGCTGTACTACGGAAGTGGATCTCTTTGGAGATGCCTTTGCAGTTCTC  
CTGGGGAGGCCCTGCAGCATCCGAAGGGGCTACCGCACCAGCTACCCCGGCCAGTGGCTGCAGCTCT  
TGATGCATGCTCAGGAAATGACCCTTTTGGCCATCTGAAGGTAGCGCAGAGGCTGCACCTGAGCTGGAC  
CTCTTTGCAATGAAGCCACCTGAGACCAGCGCTCTGTAGTTACCCCTACAGCTAGCACAGCCCTCCAG  
TTCCCGCAACTGCTCCTTCTCCTGCTCCCACGGCTGTGGCAGCCACTGCTGCCACCACCACCGCCCGC  
TGCAGCTACCACCCTGCCACCCTCTGCTGCTGCTGCCACCACCGCCGCTGCTCCTCCTGCTCTAGAT  
ATCTTTGGTGATTTGTTTATTCTGCTCCTGAAGTTGCTGCAGCACCTAAGCCAGACGGCTCCTAGCA  
TAGACCTGTTGGCAGAGATGCTTCTCCTCCCGCCACGAGGGGCTCTCCGGTGCCTGAGAGTTCTCT  
CACTGCTGACCTTTATCTGTGGACGATTTGCAGCGCCGCTCCTGCTGATCCACTGCCTCTCCTGCAAG  
GCGGAGTCTCGGGTGTATAGACCTTTTTGGGGATGCGTTTGGAAAGTGGTGTCTCTGAAACCCAGCCAG  
CACCACAGGCTGTTTCTAGTTTATCAGCATCGGCAGATCTACTAGCTGGATTTGGGGTTCTTTATGGC  
CCCTTCTACAACGCCAGTGACTCCAGCTCAGAATAACCTGCTGCAACCCAGTTTCGAGGCAGTTTTGGA  
ACGACGCCCTCGACTTCAAGCAGCAGCTTTTTGACCCATCAGGCATCTTCTGATGCCAACCATGGCAC  
CATCCGGGCAGCCTGCCCTGTCTCAATGGTCCCACCCAGTCTGCAATGGCAGCCAGCAAAGGCCCTCG  
AAGTGACCTTGACTCGTCTCTGGCCAGTTTAGTAGGCAATCTTGGGATTTCTGGTACCACATCAAAAAAG  
GGAGATCTCCAGTGAATGCTGGGGAGAAAAAGCTGACTGGTGGAGCCAACCTGGCAACCGAAAGTCACTC  
CAGCCACATGGTCAGCGGTGTTCCCCGCAAGGCACTGTTCCACCAACCAGCTCAGTCCCTCCGGTGC  
CGGGGCCCGTCCGTTGGGCAACCTGGAGCAGGATTTGGAATGCCTCCTCAGGGACAGGCATGACCATG  
ATGCTCAGCAGCCAGTCATGTTTGACAGCCCATGATGAGGCCACCCTTTGGAGCTGCAGCTGTGCCCG  
GCACACAGCTTTCTCAAGCCCTACACCTGCCACTCAGAGTCCCAAGAACTCCAGCCAAGGACCCGTT  
AGCGGATCTTAACATCAAGGATTTCTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



**Plasmid Map:**


**ACCN:** NM\_001277986

**ORF Size:** 2688 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001277986.1](#), [NP\\_001264915.1](#)

**RefSeq Size:** 4353 bp

**RefSeq ORF:** 2691 bp

**Locus ID:** 20616

**UniProt ID:** [Q61548](#)

**Cytogenetics:** 9 E3.1

**MW:** 91.8 kDa

**Gene Summary:** Adaptins are components of the adaptor complexes which link clathrin to receptors in coated vesicles. Clathrin-associated protein complexes are believed to interact with the cytoplasmic tails of membrane proteins, leading to their selection and concentration. Binding of AP180 to clathrin triskelia induces their assembly into 60-70 nm coats.[UniProtKB/Swiss-Prot Function]