

## Product datasheet for **MR231419**

### Usp15 (NM\_001301628) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Usp15 (NM_001301628) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Usp15
Synonyms:	4921514G19Rik; AI327321; E430033I05Rik; Gcap; Gcap18
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MR231419 representing NM\_001301628  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**


ATGGCGGAAGCGGAGCGCGGACCTGGACACCCAGCGTTCTGACATCGGACGCTGCTCAAACCTCGC  
 TCCGGAAAGGGGACACCTGGTATCTAGTAGATAGTCGGTGGTTCAAACAGTGGAATAAATATGTTGGCTT  
 TGACAGCTGGGACAAATACCAGATGGGAGATCAAATGTCTATCCTGGACCCATCGATAACTCTGGACTT  
 CTCAAAGATGGCGACGCGCAGTCACTTAAGGAGCACCTTATTGATGAGTTGGATTACATACTCTTGCCAA  
 CTGAGGGCTGGAATAAATTGTCACTGGTACACACTGATGGAGGGCCAGGAGCCAATAGCAAGAAAGGT  
 GGTAGAACAGGGTATGTTGTAAAGCACTGCAAAGTAGAAGTCTACCTCACAGAAGTGAAGCTCTGTGAG  
 AACGGGAACATGAACAATGTTGTAAGTTCGGAGGTTTAGCAAAGCTGACACAATAGATACGATTGAGAAGG  
 AAATAAGAAAAATCTTCAATATCCAGATGAGAAGGAGGCCAGACTGTGGAACAAGTATATGAGCAACAC  
 ATTTGAACCACTGAATAAGCCAGACAGCACCATCCAGGATGCCGGCCTATACCAAGGACAGGTGTTAGTG  
 ATAGAGCAAAAAATGAAGATGGAACGTGGCCAGGGGTCCCTCTACTCCTAATGTGAAAAACTCAAATT  
 ATGTCTCCCGTCATATACTGCTTATAAGAACTATGATTATTCAGAACCTGGAAGAAACAATGAACAGCC  
 AGGCTGTGTGGCTTAAGTAACTTGGGAAATACGTGTTTCATGAACTCAGCCATTCAGTGTTTGAGCAAC  
 ACACCTCCACTTACTGAATATTTCTCAATGATAAGTACCAAGAAGAGCTGAATTTTGACAATCCCTTAG  
 GAATGAGAGGTGAAATAGCTAAATCATACGCTGAACTCATCAAGCAGATGTGGTCTGGAAAGTTTAGCTA  
 CGTCACTCCAAGAGCATTCAAGACACAGGTAGGGCGCTTCGCACCTCAGTTCTCTGGATATCAACAGCAG  
 GACTGCCAAGAACTACTGGCTTTCTGTTGGATGGACTGCACGAGGACCTGAATAGAATTAGGAAAAAAC  
 CCTATATCCAGTTAAAAGATGCTGATGGCAGGCCAGACAAGGTGGTTGCTGAAGAAGCTGGGGAACCA  
 TTTAAAGCGAAATGATTCATCATAGTAGACATATTCATGGCCTTTTCAAGTCTACGTTAGTGTCTCT  
 GAGTGTGCTAAGATCTCGGTGACCTTTGACCCTTTCTGTTACTTGACTCTTCCATTGCCTATGAAAAAG  
 AACGTAGCTTGGAAAGTTATTTAGTTAGAATGGATCCACTTGCTAAACCTATGCAATACAAAGTTATTGT  
 GCCAAAATTGGAAATATACTAGATCTTTGTACAGCGCTGTCTGCTTTGTCGGGAGTACCTGCAGATAAG  
 ATGATAGTCACTGACATATACAATCATAGATTTTACAGGATATTTGCCGTGGATGAAAACCTGAGTAGTA  
 TTATGGAGCGGGATGATATTTACGTGTTTGAATTAACATCAACAGGGCAGAAGACACAGAGCATGTAGT  
 TATCCCTGTCTGCCTAAGAGAAAAATTCAGACACTCAAGTTACACTCACCACTGGCTCTTCACTGTTT  
 GGCCAGCCTTTTCTTATGGCTATACCACGAAACAATACTGAAGATAAACTCTATAATCTTCTACTCTTGA  
 GAATGTGCCGGTATGTCAAATGTCTACTGAAACTGAGGAAACAGATGGACACCTACGTTGCTGTGAGGA  
 CCAGAATATTAATGGGAATGGGCCAAATGGCCTACATGAAGAGGGTTACCAAGTGAGATGGAAACCGAC  
 GAGCCAGATGACGAGTCCAGCCAGGATCAAGAGCTTCCCTCAGAGAACGGAACAGTCACTGTAAGATT  
 CAGTCCGAGGAGATAACGATTCTGAAAATGGATTGTGCACTGAAGAAACATGCAAAGGTCAACTACGGG  
 ACACAAAAGCGCTTGTTTACATTCCAGTTCAACAACCTTAGGCAACAATGACATCAACTACATCAAAGAC  
 GACACCAGTCACATCCGATTTCGATGACAGGCAGCTTCGGCTGGACGAAAGATCTTTTCTCGCTTTGGACT  
 GGGACCCTGACTTGAAGAAGAGATACTTTGATGAGAATGCTGCAGAGGATTTTGA AAAACATGAAAGTGT  
 GGAATACAAACCTCCAAAAGGCCCTTTGTGAAGTTAAAAGACTGCATTGAGCTCTTACCACCAAGGAG  
 AAGCTGGGTGCTGAGGACCCCTGGTATTGTCAAATTGTAAGAGCATCAGCAAGCCACAAGAAATTAG  
 ATCTGTGGTCCCTGCCTCCAGTCTTGTGGTGCATCTCAAGAGATTCTCCTATAGTCGCTACATGAGAGA  
 CAAGCTAGATACATTAGTCGATTTTCCCATCAGTACTTGGATATGTCAGAATTCCTAATTAACCCAAAT  
 GCAGGTCTTGCCGCTATAATTTGATTGCTGTTTCCAACCACTATGGAGGGATGGGAGGAGGACATTATA  
 CTGCTTTTGCAAAAACAAGATGATGGGAAGTGGTACTATTTGATGACAGTAGCGTCTCCACTGCATC  
 TGAAGACCAAATCGTGTCAAGGCAGCGTACGTGCTCTTCTACCAGAGACAGGACACTTTCAGTGGGACT  
 GGCTTCTTCCCTCTTGACCGAGAACTAAAGGTGCTTCAGCTGCCACAGGTATCCCCTGGAAAGTGACG  
 AAGACAGCAATGATAATGACAATGACCTAGAAAATGAAAAGTGTATGCACACTAAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<a href="#">NM_001301628.1</a> , <a href="#">NP_001288557.1</a>
<b>RefSeq Size:</b>	11477 bp
<b>RefSeq ORF:</b>	2859 bp
<b>Locus ID:</b>	14479
<b>UniProt ID:</b>	<a href="#">Q8R5H1</a>
<b>Cytogenetics:</b>	10 D2
<b>MW:</b>	109.7 kDa
<b>Gene Summary:</b>	The protein encoded by this gene is a member of the large ubiquitin specific protease (Usp) family of proteins. These proteins are known to cleave ubiquitin, and contain a conserved cysteine residue (Cys box) and two conserved histidine residues (His box) that are thought to form part of the active site of the protease. This protein has been shown to cleave both the ubiquitin-proline and the ubiquitin-methionine bonds in vitro. This protein is thought to regulate many cellular processes through its deubiquitination activity, including the transforming growth factor beta (TGF-beta) pathway. Cardiac-specific overexpression of the human ortholog of this gene in mice causes enlargement of the heart that is more pronounced in the atrium than in the ventricle. This gene has two pseudogenes on chromosome 14. Alternative splicing results in multiple transcript variants that encode multiple protein isoforms.[provided by RefSeq, Aug 2014]

Product images:



Circular map for MR231419