


ACTACCACCCAAAGATCAGGATCATCACTCAGATGCTGCAGTATCACAACAAGGCCCATCTGCTCAACAT
CCCCAGCTGGAAGTGGAAAGAGGGTGATGACGCAATATGCCTTGACAGCTCAAGTTGGGTTTCATAGCC
CAGAGCTGTCTGGCTCAAGGCCTCTCCACAATGCTTGCCAATCTCTTCTATGAGGTCATTCAAAGA
TTGAGGAAGACACATGGCAGAAATACTACTTGGAAAGAGTCTCCAATGAAATGTACACAGAATATCTCTC
CAGTGCCTTCGTGGGTCTGTCTCCCTACTGTTTGTGAGCTGTGTTTGTGAAGCTTAAGCTCCTGATG
ATAGCCATTGAGTACAAGTCTGCCAACAGAGAGAGCCGAATATTAATTAACCTGGGAACCACCTTAAGA
TCCAAGAAGGTACTTTAGGATTTTTTCATCGCAAGTGATGCCAAAGAAGTTAAAAGGGCATTTTTTACTG
CAAGGCCTGCATGATGACGTACAGATCCCAAAAAGAATTAATAAATGTGGCTGCAGGCGGCTGATCTAT
TTTGAAGATGAGCAGCCGCCAACCTGTCAACAAAAAACAACGTAATGGGGCATGAGGAACCTCGC
CCAACACCTCCCCGAAGCTGATGAGGCATGACCCCTTGTTAATCCTGGCAATGATCAGATTGACAAACAT
GGACTCCAATGTGAAAAAGTACGACTCCACTGGAATGTTTCACTGGTGTGCACCAAGGAGATTGAGAAA
GTCATCTTGACTCGAAGTGAAGTGCCATGACTGTCCTGAGTGGCCATGTCGTAGTCTGCATCTTTGGG
ATGTCAGCTCAGCCCTGATTGGCTCCGGAACCTGGTGTGCCACTTCGTGCTAGCAACTTCACTATCA
TGAGCTCAAACACATTGTGTTTGTGGGCTCCATTGAGTACCTCAAGAGGGAGTGGGAAACACTGCACAAC
TTCCCGAAAGTGTCCATATTGCCTGGTACACCATTAAGTCGGGCTGATTAAGGGCTGTCAACATCAACC
TCTGTGACATGTGCGTTATCCTGTGAGCAATCAGAATAATATTGATGATACTTCGCTTCAAGGACAAAGGA
ATGCATCTTGGGCTCACTCAACATCAAATCTATGCAGTTTGTGACAGCATCGGGGTCTTGCAGGCTAAT
TCCAAGGATTCACACCTCCTGGAATGGACAGATCATCACCCGACAACAGCCAGTGCACGGGATGTTAC
GCCAGCCGTCCATCACAACCTGGGGTCAACATCCCATCATCACGGAACCTCGTGAATGATACCAATGTTCA
GTTTTTGGACCAAGACGATGACGATGACCCTGACACAGAGCTGTACCTCACACAGCCCTTTGCTTGTGGG
ACAGCATTTGCCGTGAGCGTCTGGACTCACTCATGAGCGCGACATACTTCAATGACAATATCCTCACCC
TAATACGGACCCTGGTGACAGGAGGACACACCAGAGCTCGAGGCTCTAATAGCTGAGGAGAATGCACT
TCGAGGAGGCTACAGCACTCCGACAGATTGGCCAACAGGGACCCTTCCGAGTGGCCAGTTAGCCCTG
TTAGATGGTCCCTTTGCAGACTTAGGGGATGGTGGTTGTTATGGTGATCTGTTCTGCAAAGCTCTGAAAA
CATATAATATGCTTTGTTTTGGAATTTACCGGCTGAGAGATGCCACCTCAGCACCCCGCCAGTGTAC
AAAAAGGTACGTCATACCAATCCTCCCTACGAGTTTGTGCTCGTACCAACAGACCTGATCTTCTGCCTG
ATGCAGTTTGACCACAACGCTGGCCAATCCCGGGCCAGTCTGTCTCATTCCCTCCACTCCTCACAGTCGT
CCAGTAAGAAGAGCTCCTCCGTCCACTCCATCCCGTCCACAGCAAATCGGCCGAACCGGCCAAGTCCAG
GGAGTCCCGCACAACAGAATGCAACAAGGATGACTAGAATGGGCCAAGCAGAAAAGAAATGGTTTACA
GATGAGCCGGATAATGCCTATCCAGAAACATTCAAATCAAGCCCATGAGTACCCACATGGCTAACCCAGA
TCAACCAATATAAATCCACAAGCAGCCTGATCCCTCCAATCAGAGAAGTTGAAGATGAATGT

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:



ACCN: NM_001253365

ORF Size: 3702 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001253365.1](#), [NP_001240294.1](#)

RefSeq Size: 4982 bp

RefSeq ORF: 3705 bp

Locus ID: 16531

UniProt ID: [Q08460](#)

Cytogenetics: 14 A3

MW: 137.9 kDa

Gene Summary: Potassium channel activated by both membrane depolarization or increase in cytosolic Ca(2+) that mediates export of K(+). It is also activated by the concentration of cytosolic Mg(2+). Its activation dampens the excitatory events that elevate the cytosolic Ca(2+) concentration and/or depolarize the cell membrane. It therefore contributes to repolarization of the membrane potential. Plays a key role in controlling excitability in a number of systems, such as regulation of the contraction of smooth muscle, the tuning of hair cells in the cochlea, regulation of transmitter release, and innate immunity. In smooth muscles, its activation by high level of Ca(2+), caused by ryanodine receptors in the sarcoplasmic reticulum, regulates the membrane potential. In cochlea cells, its number and kinetic properties partly determine the characteristic frequency of each hair cell and thereby helps to establish a tonotopic map. Kinetics of KCNMA1 channels are determined by alternative splicing, phosphorylation status and its combination with modulating beta subunits. Highly sensitive to both iberiotoxin (IbTx) and charybdotoxin (CTX).[UniProtKB/Swiss-Prot Function]