


Product datasheet for **MWM1001**

MYC/DDK Tagged Western Blot Molecular Weight Markers

Product data:

Product Type:	Molecular Tools
Description:	Pre-mixed molecular weight markers with 6 over-expressed cell lysates with Myc/DDK tags covering 28 to 158 KDa. Each Vial contains 50ug total protein per tube; Sufficient for 20 loads.
Formulation:	SDS Sample Buffer
Detection Antibodies:	A7, Anti-Myc tag monoclonal antibody (TA100010); DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

Product images:



Western Blot Image: 2.5ul pre-mixed molecular weight marker (MWM1001) was loaded onto the precast gel. Western blot analysis was performed with anti-DDK antibody (Cat# [TA50011-100]) at 1:2,000 dilution.



[View online »](#)