

## Product datasheet for **NP100053**

### HotStart PCR Kit

#### Product data:

Product Type:	Cloning Accessories
Applications:	Long Range PCR, High-Fidelity PCR, Fast PCR , Multiplex PCR , Specialty PCR , PCR , Routine PCR
Heat Inactivation:	No
5' - 3' Exonuclease:	No
3' - 5' Exonuclease:	Yes
Components:	HS DNA Polymerase (2 , 000U/mL) 5X HF Buffer 2.5X GC Buffer dNTPs (10mM each)
Stability:	Stable at -20°C for 1 year in the supplied storage buffer



[View online »](#)