

Product datasheet for **NP100060**

dNTPs (10 mM each)

Product data:

Product Type:	Cloning Accessories
Applications:	Whole Genome Amplification, Polymerases for DNA Manipulation, PCR, DNA Sequencing, Isothermal Amplification, Hot Start PCR, RT-PCR & cDNA Synthesis, Long Range PCR, Fast PCR, Methylated DNA Analysis, Multiplex PCR, PCR, Specialty PCR, Routine PCR.
Components:	dNTPs (10 mM each)
Stability:	Stable at -20°C for 1 year in the supplied storage buffer



[View online »](#)