

Product datasheet for **PBL21100-2**

Interferon alpha (IFN alpha) Mouse Monoclonal Antibody [Clone ID: MMHA-2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	MMHA-2
Applications:	ELISA, FN, IHC, IP, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Conjugation:	Unconjugated
Synonyms:	B Cell Interferon, IFN-alpha, IFNA, Interferon-alpha, Leukocyte Interferon, Type I Interferon



[View online »](#)