

Product datasheet for PH300163

ZNF692 (NM_017865) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	ZNF692 MS Standard C13 and N15-labeled recombinant protein (NP_060335)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC200163
Predicted MW:	57 kDa
Protein Sequence:	>RC200163 protein sequence Red =Cloning site Green =Tags(s)

MASSPAVDVSCRRREKRRQLDARRSKCRIRLGGHMEQWCLLKERLGFSLHSQLAFLLDRTYSSGCVLCA
 GPEPLPPKGLQYLVLLSHAHSRECSLVPGLRPGGQDGGLVWECSAGHTFSWGPSLSPTPSEAPKPASLP
 HTRRSWCSEATSGQELADLESEHDERTQEARLPRRVGPPPETFPPPGEEEEEDNDEDEEMLSDAS
 LWTYSSSPDDSEPDAPRLLPSPVTCTPKEGETPPAPAALSSPLAVPALSSSLSSRAPPPAEVRVQPQLS
 RTPQAAQQTEALASTGSQAQSAPTPAWDEDTAQIGPKRIRKAAKRELMPCDFPGCGRIFSNRQYLNHHKK
 YQHIHQKSFSCPEPACGKSFNFKKHLKEHMKLHSDTRDYICEFCARFRTSSNLVIHRRITGEKPLQCE
 ICGFTCRQKASLNWHQRKHAETVAALRFPCEFCGKRFEPDSVAHRKSHPALLLAPQESPSGPLEPCP
 SISAPGPLGSSEGSRPSASPQAPTLLPQQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_060335</u>
RefSeq Size:	1992
RefSeq ORF:	1557


[View online »](#)

Synonyms:	AREBP; Zfp692
Locus ID:	55657
UniProt ID:	<u>Q9BU19</u>
Cytogenetics:	1q44
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified ZNF692 protein (Cat# [TP300163]). The protein was produced from HEK293T cells transfected with ZNF692 cDNA clone (Cat# [RC200163]) using MegaTran 2.0 (Cat# [TT210002]).