


Cytogenetics: 16q22.1

Summary: This gene encodes an anti-apoptotic protein that has been shown to down-regulate the enzyme activities of caspase 2, caspase 8 and tumor protein p53. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]

Product images:



Coomassie blue staining of purified NOL3 protein (Cat# [TP300591]). The protein was produced from HEK293T cells transfected with NOL3 cDNA clone (Cat# [RC200591]) using MegaTran 2.0 (Cat# [TT210002]).