

Product datasheet for PH300656

S100A11 (NM_005620) Human Mass Spec Standard

Product data:


Product Type:	Mass Spec Standards
Description:	S100A11 MS Standard C13 and N15-labeled recombinant protein (NP_005611)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC200656
Predicted MW:	11.7 kDa
Protein Sequence:	>RC200656 protein sequence Red=Cloning site Green=Tags(s) MAKISSPTETERCIESLIAVFQKYAGKDGNYTLSKTEFLSFMNTELA AFTKNQKDPGVLD RMMKKLDTN SDGQLDFSEFLNLIIGGLAMACHDSFLKAVPSQKRT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_005611
RefSeq Size:	595
RefSeq ORF:	315
Synonyms:	HEL-S-43; MLN70; S100C
Locus ID:	6282
UniProt ID:	P31949 , V9HWH9
Cytogenetics:	1q21.3



[View online »](#)

Summary:

The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein may function in motility, invasion, and tubulin polymerization. Chromosomal rearrangements and altered expression of this gene have been implicated in tumor metastasis. [provided by RefSeq, Jul 2008]

Product images:

Coomassie blue staining of purified S100A11 protein (Cat# [TP300656]). The protein was produced from HEK293T cells transfected with S100A11 cDNA clone (Cat# [RC200656]) using MegaTran 2.0 (Cat# [TT210002]).