

## Product datasheet for PH300769

### FAM110A (NM\_001042353) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	FAM110A MS Standard C13 and N15-labeled recombinant protein (NP_001035812)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC200769
Predicted MW:	31.3 kDa
Protein Sequence:	>RC200769 protein sequence Red=Cloning site Green=Tags(s)  MPVHTLSPGAPSAPALPCRLRTRVPGYLLRGPADGGARKPSAVERLEADKAKYVKSLHVANTRQEPVQPL LSKQPLFSPETRRTVLTPSRRALPGPCRRPQLDLIDLSSLIDLCDSPVSPAEEASRTPGRAEGAGRPPAT PPRPPSTSAVRRVDVRPLPASPAPCSPGPAASSPARPPGLQRKSDLSERF SRAAADLERFFNFCG LDPEEARGLGVAHLARASSDIVSLAGPSAGPGSSEGGCSRRSSVTVEERARERVPYGVSVVERNARVIKW LYGLRQARESAAEG  TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>4</sub> ]-L-Arginine and [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>2</sub> ]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<a href="#">NP_001035812</a>
RefSeq Size:	1827
RefSeq ORF:	885
Synonyms:	bA371L19.3; C20orf55; F10
Locus ID:	83541




[View online »](#)

UniProt ID: [Q9BQ89](#)

Cytogenetics: 20p13

### Product images:



Coomassie blue staining of purified FAM110A protein (Cat# [TP300769]). The protein was produced from HEK293T cells transfected with FAM110A cDNA clone (Cat# [RC200769]) using MegaTran 2.0 (Cat# [TT210002]).