

Product datasheet for PH300951

LRRC61 (NM_023942) Human Mass Spec Standard

Product data:


Product Type:	Mass Spec Standards
Description:	LRRC61 MS Standard C13 and N15-labeled recombinant protein (NP_076431)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC200951
Predicted MW:	28 kDa
Protein Sequence:	>RC200951 protein sequence Red=Cloning site Green=Tags(s) MDPPAEKPGEAGGLQITPQLKSRTEGFSLESILLKLRGLGLADLGCLGECLGLEWLDL SGNALTHLGP LASLRQLAVLNVSNNRLTGLEPLATCENLQSLNAAGNLLATPGQLQCLAGLPCLEYLRRLRDLARLSNPL CANPSYWAAVRELLPGLKVIDGERVIGRGSEFYQLCRDL DSSLRPSSSPGPRATEAQPWVEPGYWESWPS RSSSILEEACRQFQDTLQECWDLDRQASDSLAQAEQVLSAGPTSSFVF TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_076431
RefSeq Size:	1809
RefSeq ORF:	777
Synonyms:	HSPC295
Locus ID:	65999
UniProt ID:	Q9BV99 , A0A090N7W5



[View online »](#)

Cytogenetics: 7q36.1

Product images:



Coomassie blue staining of purified LRRC61 protein (Cat# [TP300951]). The protein was produced from HEK293T cells transfected with LRRC61 cDNA clone (Cat# [RC200951]) using MegaTran 2.0 (Cat# [TT210002]).