




RefSeq Size: 9338  
RefSeq ORF: 2247  
Locus ID: 5991  
UniProt ID: [P48380](#)  
Cytogenetics: 9p24.2

**Summary:** This gene is a member of the regulatory factor X gene family, which encodes transcription factors that contain a highly-conserved winged helix DNA binding domain. The protein encoded by this gene is structurally related to regulatory factors X1, X2, X4, and X5. It is a transcriptional activator that can bind DNA as a monomer or as a heterodimer with other RFX family members. Multiple transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Aug 2013]

**Protein Families:** Transcription Factors

### Product images:



Coomassie blue staining of purified RFX3 protein (Cat# [TP301137]). The protein was produced from HEK293T cells transfected with RFX3 cDNA clone (Cat# [RC201137]) using MegaTran 2.0 (Cat# [TT210002]).