

## Product datasheet for PH302516

### UBE2H (NM\_003344) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	UBE2H MS Standard C13 and N15-labeled recombinant protein (NP_003335)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC202516
Predicted MW:	20.7 kDa
Protein Sequence:	>RC202516 protein sequence Red=Cloning site Green=Tags(s)  MSSSPSPGKRRMDTDVVKLIIESKHEVTILGGLNEFVVKFYGPQGTPEYEGVWVKVRVDLPDKYPFKSPSIGF MNKIFHPNIDEASGTVCLDVINQWTALYDLTNIFESFLPQLLAYPNPIDPLNGDAAAMYLRPEEYKQK IKEYIQKYATEEALKEQEEGTGDSSSESSMSDFSEDEAQDMEL  TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<a href="#">NP_003335</a>
RefSeq Size:	5174
RefSeq ORF:	549
Synonyms:	E2-20K; GID3; UBC8; UBCH; UBCH2
Locus ID:	7328
UniProt ID:	<a href="#">P62256</a> , <a href="#">A4D1L5</a>



[View online »](#)

**Cytogenetics:** 7q32.2

**Summary:** The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. The encoded protein sequence is 100% identical to the mouse homolog and 98% identical to the frog and zebrafish homologs. Three alternatively spliced transcript variants have been found for this gene and they encode distinct isoforms. [provided by RefSeq, Feb 2011]

**Protein Pathways:** Ubiquitin mediated proteolysis

### Product images:



Coomassie blue staining of purified UBE2H protein (Cat# [TP302516]). The protein was produced from HEK293T cells transfected with UBE2H cDNA clone (Cat# [RC202516]) using MegaTran 2.0 (Cat# [TT210002]).