

Product datasheet for PH303023

CORO1B (NM_020441) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	CORO1B MS Standard C13 and N15-labeled recombinant protein (NP_065174)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC203023
Predicted MW:	54.2 kDa
Protein Sequence:	>RC203023 protein sequence Red=Cloning site Green=Tags(s)

MSFRKVVVRQSKFRHVFGQPVKNDQCYEDIRVSRVTWDSTFCVAVNPKFLAVIVEASGGGAFLVPLSKTGR
IDKAYPTVCGHTGPVLDIDWCPHNDEVIASGSEDCTVMVWQIPENGLTSPLTEPVVLEGHTRVGIIAW
HPTARNVLLSAGCDNVVLIWNVGTAEELYRLDSLHPDLIYNVSWNHNGSLFCSACKDKSVRIIDPRRGTL
VAEREKAHEGARMRAIFLADGKVFTTGF SRMSERQLALWDPENLEEPMALQELDSSNGALLPFYDPDTS
VVVYCGKGDSSIRYFEITEEPPIHFLNTFTSKEPQRGMGSMKRGLEVSKCEIARFYKLHERKCEPIVM
TVPRKSDLFQDDLYPDAGPEAALEAEWVSGRDADPILISLREAYVPSKQRDLKISRRNVLSDSRPAMA
PGSSHLGAPASTTTAADATPSGLARAGEAGKLEEVMQELRALRALVKEQGDRICRLEEQLGRMENGDA

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_065174</u>
RefSeq Size:	1927
RefSeq ORF:	1467
Synonyms:	CORONIN-2



[View online »](#)


Locus ID: 57175

UniProt ID: [Q9BR76](#), [A0A024R5K1](#)

Cytogenetics: 11q13.2

Summary: Members of the coronin family, such as CORO1B, are WD repeat-containing actin-binding proteins that regulate cell motility (Cai et al., 2005 [PubMed 16027158]).[supplied by OMIM, Mar 2008]

Product images:



Coomassie blue staining of purified CORO1B protein (Cat# [TP303023]). The protein was produced from HEK293T cells transfected with CORO1B cDNA clone (Cat# [RC203023]) using MegaTran 2.0 (Cat# [TT210002]).