

## Product datasheet for PH303286


### SNHG32 (NM\_001040438) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	C6orf48 MS Standard C13 and N15-labeled recombinant protein (NP_001035528)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC203286
Predicted MW:	8.3 kDa
Protein Sequence:	>RC203286 protein sequence Red=Cloning site Green=Tags(s)  MERSFVWLSCLDSDSCNLTFRLLGEVESHACSPSLLNWLLTQYLPFGAGHILRTYFPVLSGVSSCHLIGG KMPEN  TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>4</sub> ]-L-Arginine and [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>2</sub> ]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<a href="#">NP_001035528</a>
RefSeq Size:	751
RefSeq ORF:	225
Synonyms:	D6S57; G8
Locus ID:	50854
Cytogenetics:	6p21.33



[View online »](#)

**Product images:**

Coomassie blue staining of purified C6orf48 protein (Cat# [TP303286]). The protein was produced from HEK293T cells transfected with C6orf48 cDNA clone (Cat# [RC203286]) using MegaTran 2.0 (Cat# [TT210002]).