

Product datasheet for PH303454

PUF60 (NM_014281) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	PUF60 MS Standard C13 and N15-labeled recombinant protein (NP_055096)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC203454
Predicted MW:	58.2 kDa
Protein Sequence:	>RC203454 protein sequence Red=Cloning site Green=Tags(s)

MATATIALQVNGQQGGGSEPA AAAAVVAAGDKWKPPQGTDSIKMENGQSTAAKLG LPLTPEQQEALQKA
KKYAMEQSIKSVLVKQ TIAHQQQLTNLQMAAQRQRALAIMCRVYVGSIIYELGEDTIRQAFAPFGPIKS
IDMSWDSVTMKHKGF AFVEYEVPEAAQLALEQMNSVMLGGRNIKVGRPSNIGQAQPIIDQLAEEARAFNR
IYVASVHQDLSDDDIKSVFEAFGKIKSCTLARDPTTGKHKGYGFIEYEKAQSSQDAVSSMNLFDLGGQYL
RVGKAVTPPMP L LTPATPGGLPPAAAVAAAAATAKITAQEAVAGAAVLGTLGTPGLVSPALTLAQPLGTL
PQAVMAAQAPGVITGVTPARPPIPVTIPSVGVVNPILASPPTLGLLEPKKEKEEEELFPESERPEMLSEQ
EHMSISGSSARHMVMQKLLRKQESTVMVLRNMVDPKDIDDDLEGEVTEECGKFGAVNRV I IYQEKQGEE
DAEIIVKIFVEFSIASETHKAIQALNGRWFAGRKVVAEVYDQERFDNSDLSA

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_055096
RefSeq Size:	1897
RefSeq ORF:	1626



[View online »](#)

Synonyms: FIR; RoBPI; SIAHBP1; VRJS

Locus ID: 22827


UniProt ID: [Q9UHX1](#)

Cytogenetics: 8q24.3

Summary: This gene encodes a nucleic acid-binding protein that plays a role in a variety of nuclear processes, including pre-mRNA splicing and transcriptional regulation. The encoded protein forms a complex with the far upstream DNA element (FUSE) and FUSE-binding protein at the myelocytomatosis oncogene (MYC) promoter. This complex represses MYC transcription through the core-TFIID basal transcription factor. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2012]

Protein Pathways: Spliceosome

Product images:



Coomassie blue staining of purified PUF60 protein (Cat# [TP303454]). The protein was produced from HEK293T cells transfected with PUF60 cDNA clone (Cat# [RC203454]) using MegaTran 2.0 (Cat# [TT210002]).