


Product datasheet for PH303572

C6orf136 (NM_145029) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	C6orf136 MS Standard C13 and N15-labeled recombinant protein (NP_659466)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC203572
Predicted MW:	20.5 kDa
Protein Sequence:	>RC203572 protein sequence Red=Cloning site Green=Tags(s) MEEHLSVMYERLRQELPKLFLQSHDYSLSLDVEFINEILNIRTKGRTWYILSLTLCRFLAWNYFAHLRL EVLQLTRHPENWTLQARWRLVGLPVHLLFLRFYKRDKDEHYRTYDAYSTFYLNSSGLICRHRDLKLMPSH SPPTPVKKLLVGALVALGLSEPEPDNLCSKP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_659466
RefSeq Size:	1039
RefSeq ORF:	516
Locus ID:	221545
UniProt ID:	Q5SQH8
Cytogenetics:	6p21.33


[View online »](#)

Product images:


Coomassie blue staining of purified C6orf136 protein (Cat# [TP303572]). The protein was produced from HEK293T cells transfected with C6orf136 cDNA clone (Cat# [RC203572]) using MegaTran 2.0 (Cat# [TT210002]).