

## Product datasheet for PH303617

### LRRC46 (NM\_033413) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	LRRC46 MS Standard C13 and N15-labeled recombinant protein (NP_219481)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC203617
Predicted MW:	35.3 kDa
Protein Sequence:	>RC203617 protein sequence Red=Cloning site Green=Tags(s)

MSGGKSAQGPEEGVCITEALITKRNLTFPEDGELSEKMFHTLDELQTVRLDREGITTIRNLEGLQNLHS  
LYLQGNKIQQIENLACIPSLRFLSLAGNQIRQVENLLDLPLCLQFLDLSENLIETLKLDEFQPSLLILNLS  
GNSCTNQDGYRELVTEALPLLLDLGQPVVERWISDEEDEASSDEEFPELSPGFCSEKRGFLKELEQELSR  
HREHRQQTALTEHLLRMEMQPTLTDLPLPGVPMAGDSSPSATPAQGEETVPEAVSSPQASSPTKKPCSL  
IPRGHQSSFWGRKGARAATAPKASVAEAPSTTKTTAKRSKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>4</sub> ]-L-Arginine and [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>2</sub> ]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u><a href="#">NP_219481</a></u>
RefSeq Size:	1938
RefSeq ORF:	963
Locus ID:	90506
UniProt ID:	<u><a href="#">Q96FV0</a></u>



[View online »](#)

Cytogenetics: 17q21.32

### Product images:



Coomassie blue staining of purified LRRC46 protein (Cat# [TP303617]). The protein was produced from HEK293T cells transfected with LRRC46 cDNA clone (Cat# [RC203617]) using MegaTran 2.0 (Cat# [TT210002]).