

Product datasheet for PH303808

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PYK2 (PTK2B) (NM 004103) Human Mass Spec Standard

Product data:

Product Type: Mass Spec Standards

Description: PTK2B MS Standard C13 and N15-labeled recombinant protein (NP_004094)

Species: Human **HEK293 Expression Host: Expression cDNA Clone**

or AA Sequence:

RC203808

Predicted MW: 115.9 kDa

>RC203808 protein sequence **Protein Sequence:**

Red=Cloning site Green=Tags(s)

MSGVSEPLSRVKLGTLRRPEGPAEPMVVVPVDVEKEDVRILKVCFYSNSFNPGKNFKLVKCTVQTEIREI ITSILLSGRIGPNIRLAECYGLRLKHMKSDEIHWLHPQMTVGEVQDKYECLHVEAEWRYDLQIRYLPEDF MESLKEDRTTLLYFYQQLRNDYMQRYASKVSEGMALQLGCLELRRFFKDMPHNALDKKSNFELLEKEVGL DLFFPKQMQENLKPKQFRKMIQQTFQQYASLREEECVMKFFNTLAGFANIDQETYRCELIQGWNITVDLV IGPKGIRQLTSQDAKPTCLAEFKQIRSIRCLPLEEGQAVLQLGIEGAPQALSIKTSSLAEAENMADLIDG YCRLQGEHQGSLIIHPRKDGEKRNSLPQIPMLNLEARRSHLSESCSIESDIYAEIPDETLRRPGGPQYGI AREDVVLNRILGEGFFGEVYEGVYTNHKGEKINVAVKTCKKDCTLDNKEKFMSEAVIMKNLDHPHIVKLI GIIEEEPTWIIMELYPYGELGHYLERNKNSLKVLTLVLYSLQICKAMAYLESINCVHRDIAVRNILVASP ECVKLGDFGLSRYIEDEDYYKASVTRLPIKWMSPESINFRRFTTASDVWMFAVCMWEILSFGKQPFFWLE NKDVIGVLEKGDRLPKPDLCPPVLYTLMTRCWDYDPSDRPRFTELVCSLSDVYQMEKDIAMEQERNARYR TPKILEPTAFQEPPPKPSRPKYRPPPQTNLLAPKLQFQVPEGLCASSPTLTSPMEYPSPVNSLHTPPLHR HNVFKRHSMREEDFIQPSSREEAQQLWEAEKVKMRQILDKQQKQMVEDYQWLRQEEKSLDPMVYMNDKSP LTPEKEVGYLEFTGPPOKPPRLGAQSIQPTANLDRTDDLVYLNVMELVRAVLELKNELCQLPPEGYVVVV KNVGLTLRKLIGSVDDLLPSLPSSSRTEIEGTQKLLNKDLAELINKMRLAQQNAVTSLSEECKRQMLTAS HTLAVDAKNLLDAVDQAKVLANLAHPPAE

TRTRPLEQKLISEEDLAANDILDYKDDDDK**V**

C-Myc/DDK Tag:

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Concentration: >0.05 µg/µL as determined by microplate BCA method

Labeling Method: Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine

25 mM Tris-HCl, 100 mM glycine, pH 7.3 **Buffer:**

Storage: Store at -80°C. Avoid repeated freeze-thaw cycles.





Summary:

PYK2 (PTK2B) (NM_004103) Human Mass Spec Standard - PH303808

Stability: Stable for 3 months from receipt of products under proper storage and handling conditions.

RefSeq: NP 004094

RefSeq Size: 4193 RefSeq ORF: 3027

Synonyms: CADTK; CAKB; FADK2; FAK2; PKB; PTK; PYK2; RAFTK

Locus ID: 2185 UniProt ID: Q14289 Cytogenetics: 8p21.2

> This gene encodes a cytoplasmic protein tyrosine kinase which is involved in calcium-induced regulation of ion channels and activation of the map kinase signaling pathway. The encoded protein may represent an important signaling intermediate between neuropeptide-activated receptors or neurotransmitters that increase calcium flux and the downstream signals that regulate neuronal activity. The encoded protein undergoes rapid tyrosine phosphorylation and activation in response to increases in the intracellular calcium concentration, nicotinic acetylcholine receptor activation, membrane depolarization, or protein kinase C activation. This protein has been shown to bind CRK-associated substrate, nephrocystin, GTPase regulator associated with FAK, and the SH2 domain of GRB2. The encoded protein is a member of the FAK subfamily of protein tyrosine kinases but lacks significant sequence similarity to kinases from other subfamilies. Four transcript variants encoding two different


isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Protein Kinase

Protein Pathways: Calcium signaling pathway, Chemokine signaling pathway, GnRH signaling pathway,

Leukocyte transendothelial migration, Natural killer cell mediated cytotoxicity

Product images:



Coomassie blue staining of purified PTK2B protein (Cat# [TP303808]). The protein was produced from HEK293T cells transfected with PTK2B cDNA clone (Cat# [RC203808]) using MegaTran 2.0 (Cat# [TT210002]).