

## **Product datasheet for PH303819**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CD14 (NM 000591) Human Mass Spec Standard

**Product data:** 

**Product Type:** Mass Spec Standards

**Description:** CD14 MS Standard C13 and N15-labeled recombinant protein (NP\_000582)

Species:HumanExpression Host:HEK293

Expression cDNA Clone

or AA Sequence:

RC203819

Predicted MW: 40.08 kDa

**Protein Sequence:** >RC203819 representing NM\_000591

Red=Cloning site Green=Tags(s)

MERASCLLLLLPLVHVSATTPEPCELDDEDFRCVCNFSEPQPDWSEAFQCVSAVEVEIHAGGLNLEPFL KRVDADADPRQYADTVKALRVRRLTVGAAQVPAQLLVGALRVLAYSRLKELTLEDLKITGTMPPLPLEAT GLALSSLRLRNVSWATGRSWLAELQQWLKPGLKVLSIAQAHSPAFSCEQVRAFPALTSLDLSDNPGLGER GLMAALCPHKFPAIQNLALRNTGMETPTGVCAALAAAGVQPHSLDLSHNSLRATVNPSAPRCMWSSALNS LNLSFAGLEQVPKGLPAKLRVLDLSCNRLNRAPQPDELPEVDNLTLDGNPFLVPGTALPHEGSMNSGVVP

ACARSTLSVGVSGTLVLLQGARGFA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Concentration:** >0.05 μg/μL as determined by microplate BCA method

Labeling Method: Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3

**Storage:** Store at -80°C. Avoid repeated freeze-thaw cycles.

**Stability:** Stable for 3 months from receipt of products under proper storage and handling conditions.

**RefSeq:** NP 000582

 RefSeq Size:
 1623

 RefSeq ORF:
 1125

 Locus ID:
 929

 UniProt ID:
 P08571





Cytogenetics: 5q31.3

**Summary:** The protein encoded by this gene is a surface antigen that is preferentially expressed on

> monocytes/macrophages. It cooperates with other proteins to mediate the innate immune response to bacterial lipopolysaccharide, and to viruses. This gene has been identified as a target candidate in the treatment of SARS-CoV-2-infected patients to potentially lessen or inhibit a severe inflammatory response. Alternative splicing results in multiple transcript

variants encoding the same protein. [provided by RefSeq, Aug 2020]


**Protein Families:** Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS,

Transmembrane

**Protein Pathways:** Hematopoietic cell lineage, MAPK signaling pathway, Pathogenic Escherichia coli infection,

Regulation of actin cytoskeleton, Toll-like receptor signaling pathway

## **Product images:**



Coomassie blue staining of purified CD14 protein (Cat# [TP303819]). The protein was produced from HEK293T cells transfected with CD14 cDNA clone (Cat# [RC203819]) using MegaTran 2.0

(Cat# [TT210002]).