

Product datasheet for PH303850

C9orf30 (MSANTD3) (NM_080655) Human Mass Spec Standard

Product data:


Product Type:	Mass Spec Standards
Description:	C9orf30 MS Standard C13 and N15-labeled recombinant protein (NP_542386)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC203850
Predicted MW:	32.4 kDa
Protein Sequence:	>RC203850 protein sequence Red=Cloning site Green=Tags(s) MQNNEIIPAKYFSELEKSILLALVEKYKYVLECKSDARTIALKQRTWQALAHEYNSQPSVSLRDFKQL KKCWENIKARTKKIMAHERRKVKRSVSPLLSTHVLGKEKIASMLPEQLYFLQSPPEEEPEYHPDASAQE SFAVSNRELCDDEKEFIHFPVCEGTSQPEPSCSAVRITANKNYRSKTSQEGALKKMHEEEHHQMSILQL QLIQMNEVHVAKIQQIERECEMAEEHRIKMEVLNKKMYWERKLQFTKEWPVSSFNRFPNPS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_542386
RefSeq Size:	1826
RefSeq ORF:	828
Synonyms:	C9orf30; L8
Locus ID:	91283
UniProt ID:	Q96H12 , A0A024R171



[View online »](#)

Cytogenetics: 9q31.1

Product images:



Coomassie blue staining of purified MSANTD3 protein (Cat# [TP303850]). The protein was produced from HEK293T cells transfected with MSANTD3 cDNA clone (Cat# [RC203850]) using MegaTran 2.0 (Cat# [TT210002]).