

## Product datasheet for PH304931

### LENG1 (NM\_024316) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	LENG1 MS Standard C13 and N15-labeled recombinant protein (NP_077292)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC204931
Predicted MW:	30.5 kDa
Protein Sequence:	>RC204931 protein sequence Red=Cloning site Green=Tags(s)


MNILPKKSWHVRNKDNVARVRRDEAQAREEEKERERRVLLAQQEARTFLRKKARHQNSLPELEAAEAGA  
PGSGPVDFRELLLEEGKGVIRGNKEYEEERQEKERQEKALGILTYLGQSAAEAQTQPPWYQLPPGRGGP  
PPGPAPDEKIKSRLDPLREMQKHLGKKRQHGGDEGSRSRKEKEGSEKQRPKEPPSLDQLRAERLRREAAE  
RSRAEALLARVQGRALQEGQPEEDETDDRRRRYNSQFNPQLARRPRQQDPHLTH

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>4</sub> ]-L-Arginine and [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>2</sub> ]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<a href="#">NP_077292</a>
RefSeq Size:	941
RefSeq ORF:	792
Locus ID:	79165
UniProt ID:	<a href="#">Q96BZ8</a>
Cytogenetics:	19q13.42



[View online »](#)

**Product images:**

Coomassie blue staining of purified LENG1 protein (Cat# [TP304931]). The protein was produced from HEK293T cells transfected with LENG1 cDNA clone (Cat# [RC204931]) using MegaTran 2.0 (Cat# [TT210002]).