

Product datasheet for PH305101

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MPST (NM_021126) Human Mass Spec Standard

Product data:

Product Type: Mass Spec Standards

Description: MPST MS Standard C13 and N15-labeled recombinant protein (NP_066949)

Species:HumanExpression Host:HEK293

Expression cDNA Clone

RC205101

or AA Sequence: Predicted MW:

33.2 kDa

Protein Sequence: >RC205101 protein sequence

Red=Cloning site Green=Tags(s)

MASPQLCRALVSAQWVAEALRAPRAGQPLQLLDASWYLPKLGRDARREFEERHIPGAAFFDIDQCSDRTS PYDHMLPGAEHFAEYAGRLGVGAATHVVIYDASDQGLYSAPRVWWMFRAFGHHAVSLLDGGLRHWLRQNL PLSSGKSQPAPAEFRAQLDPAFIKTYEDIKENLESRRFQVVDSRATGRFRGTEPEPRDGIEPGHIPGTVN IPFTDFLSQEGLEKSPEEIRHLFQEKKVDLSKPLVATCGSGVTACHVALGAYLCGKPDVPIYDGSWVEWY

MRARPEDVISEGRGKTH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Concentration: >0.05 μg/μL as determined by microplate BCA method

Labeling Method: Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3

Storage: Store at -80°C. Avoid repeated freeze-thaw cycles.

Stability: Stable for 3 months from receipt of products under proper storage and handling conditions.

RefSeq: NP 066949

RefSeq Size: 1361 RefSeq ORF: 891

Synonyms: MST; TST2; TUM1

Locus ID: 4357



UniProt ID: P25325

Cytogenetics: 22q12.3

Summary: This protein encoded by this gene catalyzes the transfer of a sulfur ion from 3-

mercaptopyruvate to cyanide or other thiol compounds. It may be involved in cysteine degradation and cyanide detoxification. There is confusion in literature between this protein

(mercaptopyruvate sulfurtransferase, MPST), which appears to be cytoplasmic, and

thiosulfate sulfurtransferase (rhodanese, TST, GeneID:7263), which is a mitochondrial protein. Deficiency in MPST activity has been implicated in a rare inheritable disorder known as


mercaptolactate-cysteine disulfiduria (MCDU). Alternatively spliced transcript variants encoding same or different isoforms have been identified for this gene. [provided by RefSeq,

Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Cysteine and methionine metabolism, Metabolic pathways

Product images:



Coomassie blue staining of purified MPST protein (Cat# [TP305101]). The protein was produced from HEK293T cells transfected with MPST cDNA clone (Cat# [RC205101]) using MegaTran 2.0 (Cat# [TT210002]).