

Product datasheet for PH305215

C20orf141 (NM_080739) Human Mass Spec Standard

Product data:


Product Type:	Mass Spec Standards
Description:	C20orf141 MS Standard C13 and N15-labeled recombinant protein (NP_542777)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC205215
Predicted MW:	17.4 kDa
Protein Sequence:	>RC205215 protein sequence Red=Cloning site Green=Tags(s) MTRLCLPRPEAREDPIPVPPRGLGAGEGSGSPVPPVSTWGPSWAQLLDSVLWL GALGLTIQAVFSTTGP ALLLLLVSFLTFDLLHRPAGHTLPQRKLLTRGQSQGAGEGPGQEALLLQMGTVSGQLSLQDALLLLMG LGPLL RACGMPLTLLGLAFCLHPWA TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_542777
RefSeq Size:	741
RefSeq ORF:	495
Synonyms:	dj860F19.4
Locus ID:	128653
UniProt ID:	Q9NUB4



[View online »](#)

Cytogenetics: 20p13

Product images:



Coomassie blue staining of purified C20orf141 protein (Cat# [TP305215]). The protein was produced from HEK293T cells transfected with C20orf141 cDNA clone (Cat# [RC205215]) using MegaTran 2.0 (Cat# [TT210002]).