

Product datasheet for PH305447

ATP6V1B2 (NM_001693) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	ATP6V1B2 MS Standard C13 and N15-labeled recombinant protein (NP_001684)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC205447
Predicted MW:	56.4 kDa
Protein Sequence:	>RC205447 protein sequence Red=Cloning site Green=Tags(s)

MALRAMRGIVNGAAPELVPVTGGPAVGAQEALAVSRNYLSQPRLTYKTVSGVNGPLVILDHVKFPRYAE
IVHLTLPDGTKRSGQVLEVSGSKAVVQVFEGTSGIDAKKTSCEFTGDILRTPVSEDMLGRVFNKSGKPID
RGPVVLAEFDLIMGQPINPQCRIYPEEMIRTGISAIDGMNSIARGQKIPFSAAGLPHNEIAAQICRQA
GLVKKSKDVVDYSEENFAIVFAAMGVNMETARFFKSDFEENGSMNDVCLFLNLANOPTIERIITPRLALT
TAEFLAYQCEKHVLVILTDMSSYAEALREVSAAAREEVPGRRGFPGYMYTDLATIYERAGRVGGRNGSITQ
IPILTMPNDDITHPIPDLTGYITEGQIYVDRQLHNRQIYPPINVLPSLSRLMKSATGEGMTRKDHADVSN
QLYACYAIGKDVQAMKAVVGEEALTSDDLLEYLFLQKFERNFIAQGPYENRTVFETLDIGWQLLRIFPKE
MLKRIPQSTLSEFYPRDSAKH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001684
RefSeq Size:	3054
RefSeq ORF:	1533



[View online »](#)

Synonyms: ATP6B1B2; ATP6B2; DOOD; HO57; VATB; Vma2; VPP3; ZLS2

Locus ID: 526

UniProt ID: [P21281](#), [A0A140VK65](#)


Cytogenetics: 8p21.3

Summary: This gene encodes a component of vacuolar ATPase (V-ATPase), a multisubunit enzyme that mediates acidification of eukaryotic intracellular organelles. V-ATPase dependent organelle acidification is necessary for such intracellular processes as protein sorting, zymogen activation, receptor-mediated endocytosis, and synaptic vesicle proton gradient generation. V-ATPase is composed of a cytosolic V1 domain and a transmembrane V0 domain. The V1 domain consists of three A, three B, and two G subunits, as well as a C, D, E, F, and H subunit. The V1 domain contains the ATP catalytic site. The protein encoded by this gene is one of two V1 domain B subunit isoforms and is the only B isoform highly expressed in osteoclasts. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Epithelial cell signaling in Helicobacter pylori infection, Metabolic pathways, Oxidative phosphorylation, Vibrio cholerae infection

Product images:



Coomassie blue staining of purified ATP6V1B2 protein (Cat# [TP305447]). The protein was produced from HEK293T cells transfected with ATP6V1B2 cDNA clone (Cat# [RC205447]) using MegaTran 2.0 (Cat# [TT210002]).