

## Product datasheet for PH306078

### ZCCHC13 (NM\_203303) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	ZCCHC13 MS Standard C13 and N15-labeled recombinant protein (NP_976048)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC206078
Predicted MW:	18 kDa
Protein Sequence:	>RC206078 protein sequence Red=Cloning site Green=Tags(s)  MSSKDFFACGHSQHWARGCPRGGAGRRGGGHGRGSQCGSTTLSYTCYCCGESGRNAKNCVLLGNICYNC GRSGHIAKDCKDPKRERQHCYTCGRLGHLARDQKEQKCYSCGKLGHIQKDCAQVKCYRCGEIGHVA INCSKARPGQLLPLRQIPTSSQGMSQ  TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<a href="#">NP_976048</a>
RefSeq Size:	859
RefSeq ORF:	498
Synonyms:	CNBP2; ZNF9L
Locus ID:	389874
UniProt ID:	<a href="#">Q8WW36</a>




[View online »](#)

**Cytogenetics:** Xq13.2

**Summary:** This gene appears to represent an intronless retrocopy of a related multi-exon gene located on chromosome 3. However, the CDS of this intronless gene remains relatively intact, it is conserved in other mammalian species, it is known to be transcribed, and it is therefore thought to encode a functional protein. The encoded protein contains six CCHC-type zinc fingers, and is thus thought to function as a transcription factor. [provided by RefSeq, May 2010]

**Product images:**



Coomassie blue staining of purified ZCCHC13 protein (Cat# [TP306078]). The protein was produced from HEK293T cells transfected with ZCCHC13 cDNA clone (Cat# [RC206078]) using MegaTran 2.0 (Cat# [TT210002]).