

## Product datasheet for PH306241


### SH3BGRL2 (NM\_031469) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	SH3BGRL2 MS Standard C13 and N15-labeled recombinant protein (NP_113657)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC206241
Predicted MW:	12.1 kDa
Protein Sequence:	>RC206241 representing NM_031469 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MVirVFIASSSGFVAIKKKQDVVRFLEANKIEFEEVDITMSEEQRQWMYKNVPPEKKPTQGNLPPQIF NGDRYCGDYDSFFESKESNTVFSFLGLKPRLASKAEP  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<a href="#">NP_113657</a>
RefSeq Size:	4676
RefSeq ORF:	321
Locus ID:	83699
UniProt ID:	<a href="#">Q9UJC5</a>
Cytogenetics:	6q14.1



[View online »](#)

**Product images:**

Coomassie blue staining of purified SH3BGRL2 protein (Cat# [TP306241]). The protein was produced from HEK293T cells transfected with SH3BGRL2 cDNA clone (Cat# [RC206241]) using MegaTran 2.0 (Cat# [TT210002]).