

Product datasheet for PH306684

C1orf87 (NM_152377) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	C1orf87 MS Standard C13 and N15-labeled recombinant protein (NP_689590)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC206684
Predicted MW:	62 kDa
Protein Sequence:	>RC206684 protein sequence Red=Cloning site Green=Tags(s)

MSSAWKTPRGSDAMPEIMVKIIGSKHFQYLVEKPKIKENDSLKTETQTMHQKPMTDNARQMSRDTVPVPI
FTDQQTDDNPDDVKEKKHPENNQKSENNQKLLTGANSSRFLDGNIPSQANVHCSSVPTGDQSLSYVHGIP
RRKLRDWSLEQMVRGSSDQPEDIQSPSGTTNEDAFLLALVRRELKSRPLSSNLEKLQKELKILDPISS
GFLLQSQLSRLFLKHEVPLQLPTVKILCQRFSKRGSPPEMNYEKLLWFLNSAASDYPQQNKAAADLRKTE
SHGTHSQSTPPQHSSSQPEVNRSLLEILKMALRTTNGRLNIDNLNLSFRKEDRSFSGCLPLPKVRAICGK
HGLYLTLALLELLNHQDLGYQNEIKWQNFVEMLTRASSDLLSDLPTGKNEKKAPAPMEPEVPEMSQSK
TEHMKTPPEELQPESSPAETSACKDPLKPLKIRPVVSQPFVNPVAVKNKAAECETWIDFRFKLENALYLCDL
SNTGVLEKERARRLIHNYNLIYNLSLSPQKIDQALRRFRSGENMLLEPALRYLKEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_689590
RefSeq Size:	2126
RefSeq ORF:	1638



[View online »](#)

Synonyms: CREF
Locus ID: 127795
UniProt ID: [Q8N0U7](#)
Cytogenetics: 1p32.1

Product images:



Coomassie blue staining of purified C1orf87 protein (Cat# [TP306684]). The protein was produced from HEK293T cells transfected with C1orf87 cDNA clone (Cat# [RC206684]) using MegaTran 2.0 (Cat# [TT210002]).