

Product datasheet for PH306880

CXorf38 (NM_144970) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	CXorf38 MS Standard C13 and N15-labeled recombinant protein (NP_659407)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC206880
Predicted MW:	36.7 kDa
Protein Sequence:	>RC206880 protein sequence Red =Cloning site Green =Tags(s)

MVLSELAARLNCAEYKNWVKAGHCLLLLRSLQGFVGREVLSFHRGLAAPGLGPRAVCRGGSRCSPRA
 RQFQPQCQVCAEWKREILRHVNRNGDVHWGNCRPGRWPVDAWEVAKAFMPRGLADKQGPEECDVALLS
 LINSDFVVDKRVTEVIKCRNEIMHSSEMKVSSTWLRDFQMKIQNFLNEFKNIPEIVAVYSRIEQLLT
 SDWAVHIPEEDQRDGCECEMGTYLSESQVNEIEMQLLKEKLQEIYLQAEQEVLPEELSNRLEVVKFLR
 NNEDLRNGLTEDMQKLDLCLHQLDSQEPGRQTPDRKA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_659407</u>
RefSeq Size:	4260
RefSeq ORF:	957
Locus ID:	159013
UniProt ID:	<u>Q8TB03</u>


[View online »](#)

Cytogenetics: Xp11.4

Product images:



Coomassie blue staining of purified CXorf38 protein (Cat# [TP306880]). The protein was produced from HEK293T cells transfected with CXorf38 cDNA clone (Cat# [RC206880]) using MegaTran 2.0 (Cat# [TT210002]).