

Product datasheet for PH306897

Ubiquilin (UBQLN1) (NM_013438) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	UBQLN1 MS Standard C13 and N15-labeled recombinant protein (NP_038466)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC206897
Predicted MW:	62.6 kDa
Protein Sequence:	>RC206897 protein sequence Red=Cloning site Green=Tags(s)

MAESGESGGPPGSQDSAAGAEGAGTPAAAAAEPKIMKVTVKTPKEKEEFVAVPENSSVQQFKEEISKRFK
SHTDQLVLIFAGKILKDQDTLSQHGIDGLTVHLVIKTQNRPDHSAQQTNTAGSNVTTSSTPNSNSTSG
SATSNPFLGGLGGLAGLSSLGLNTNFSELQSQMQRQLLSNPMMVQIMENPFVQSMLSNHDLMRQLIM
ANPQMQLIQRNPEISHMLNNDIMRQTLELARNPAMMQEMMRNQDRALSNLESIPGGYNALRRMYTDIQ
EPMLSAAQEQFGGNPFASLVNNTSSGEGSQPSRTENRDPLPNPWAPQTSQSSASSTASTVGGTTGSTA
SGTSGQSTTAPNLVPGVGASMFNTPGMQSLQLQITENPQLMQNMLSAPYMRSMQSLSQNPDLAAQMMLN
NPLFAGNPQLQEQMRQLPTFLQQMQNPDTLSAMSNPRAMQALLQIQQLQLATEAPGLIPGFTPLGLGA
LGSTGGSSGTNGSNATPSENTSPTAGTTEPGHQFIIQQMLQALAGVNPQLQNPEVRFQQLEQLSAMGFL
NREANLQALIATGGDINAAIERLLGSQPS

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_038466
RefSeq Size:	4177
RefSeq ORF:	1767



[View online »](#)

Synonyms: DA41; DSK2; PLIC-1; UBQN; XDRP1

Locus ID: 29979


UniProt ID: [Q9UMX0](#)

Cytogenetics: 9q21.2-q21.3

Summary: This gene encodes an ubiquitin-like protein (ubiquilin) that shares a high degree of similarity with related products in yeast, rat and frog. Ubiquilins contain an N-terminal ubiquitin-like domain and a C-terminal ubiquitin-associated domain. They physically associate with both proteasomes and ubiquitin ligases, and thus are thought to functionally link the ubiquitination machinery to the proteasome to affect in vivo protein degradation. This ubiquilin has also been shown to modulate accumulation of presenilin proteins, and it is found in lesions associated with Alzheimer's and Parkinson's disease. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified UBQLN1 protein (Cat# [TP306897]). The protein was produced from HEK293T cells transfected with UBQLN1 cDNA clone (Cat# [RC206897]) using MegaTran 2.0 (Cat# [TT210002]).