

Product datasheet for PH307057

C13orf26 (TEX26) (NM_152325) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	C13orf26 MS Standard C13 and N15-labeled recombinant protein (NP_689538)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC207057
Predicted MW:	33.6 kDa
Protein Sequence:	>RC207057 protein sequence Red =Cloning site Green =Tags(s)

MEQPGPRAPDPSLCHHNLQPTDDPNWDSYATTMRTAFTPKTGAVPALIRQNGIRRLGYTYSLSDPILNQT
 QYSDEYTWKSHSKEDLIKTETSRGIKSHKSHLNEDIFLWTLPHCQQTGTLKNCLPWKIPASMKEVKNKALS
 NQFISLTKRDFVDRSKAQKIKKSSHLLEWKKLLPQPPDTEFRRNYQIPAKIPELQDFSFKYGCYSSLPV
 ASQGLVPSVLHSLRNQEHTKKQTTYQSDYDKTYPDFLMLLSFTSSQVKEYLQSLSYKDRQIIDRFIRT
 HCDTNKKKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_689538</u>
RefSeq Size:	1113
RefSeq ORF:	867
Synonyms:	C13orf26
Locus ID:	122046


[View online »](#)

UniProt ID: Q8N6G2

Cytogenetics: 13q12.3

Product images:



Coomassie blue staining of purified TEX26 protein (Cat# [TP307057]). The protein was produced from HEK293T cells transfected with TEX26 cDNA clone (Cat# [RC207057]) using MegaTran 2.0 (Cat# [TT210002]).