

Product datasheet for PH308472

PLAC1L (OOSP2) (NM_173801) Human Mass Spec Standard

Product data:


Product Type:	Mass Spec Standards
Description:	PLAC1L MS Standard C13 and N15-labeled recombinant protein (NP_776162)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC208472
Predicted MW:	18 kDa
Protein Sequence:	>RC208472 protein sequence Red=Cloning site Green=Tags(s) MALEVLMLLAVLIWTGAENLHVKISCSLDWLMVSVIPVAESRNLYIFADELHLGMGCPANRIHTYVYEFI YLVRCGIRTRVYSEETLLFQTELYFTPRNIDHDPQEIHLECSTSRKSVWLTPVSTENEIKLDPSPIAD FQTTAEELGLLSSPNLL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_776162</u>
RefSeq Size:	1646
RefSeq ORF:	474
Synonyms:	OOSP2A; PLAC1L; TMEM122
Locus ID:	219990
UniProt ID:	<u>Q86WS3</u> , <u>A0A140VIR6</u>


[View online »](#)

Cytogenetics: 11q12.1

Protein Families: Secreted Protein

Product images:



Coomassie blue staining of purified OOSP2 protein (Cat# [TP308472]). The protein was produced from HEK293T cells transfected with OOSP2 cDNA clone (Cat# [RC208472]) using MegaTran 2.0 (Cat# [TT210002]).