

Product datasheet for PH308488

C6orf62 (NM_030939) Human Mass Spec Standard

Product data:


Product Type:	Mass Spec Standards
Description:	C6orf62 MS Standard C13 and N15-labeled recombinant protein (NP_112201)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC208488
Predicted MW:	27.1 kDa
Protein Sequence:	>RC208488 protein sequence Red=Cloning site Green=Tags(s) MGDPNSRKKQALNRLRAQLRKKKESLADQDFDKMYIAFVFKEKKKKSALFEVSEVIPVMTNNYEENILKG VRDSSYSLESSLELLQKDVVQLHAPRYQSMRRDVIQCTQEMDFILWPRNDIEKIVCLLFSRWKESDEPFR PVQAKFEFHGDYEQFLHVLSRKDKTGIVVNNPNQSVFLFIDRQHLQTPKNKATIFKLCSICLYLPQEQ LTHWAVGTIEDHLRPYMPE TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_112201</u>
RefSeq Size:	2498
RefSeq ORF:	687
Synonyms:	dj30M3.2; Nbla00237; XTP12
Locus ID:	81688
UniProt ID:	<u>Q9GZU0, A0A024R026</u>



[View online »](#)

Cytogenetics: 6p22.3

Product images:



Coomassie blue staining of purified C6orf62 protein (Cat# [TP308488]). The protein was produced from HEK293T cells transfected with C6orf62 cDNA clone (Cat# [RC208488]) using MegaTran 2.0 (Cat# [TT210002]).