

Product datasheet for PH308543

ZNF26 (NM_019591) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	ZNF26 MS Standard C13 and N15-labeled recombinant protein (NP_062537)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC208543
Predicted MW:	61.3 kDa
Protein Sequence:	>RC208543 protein sequence Red=Cloning site Green=Tags(s)

MATSFRTASCWGLLSFKDISMEFTWDEWQLLDSTQKYLVRDVIENYHNLISVGYHGTPDLIFKLEQGE
DPWIINAKISRQSCPDGWEWYQNNQDELESIERSYACSVLGRNLKTHDSSRQLYNTRGKSLTQNSA
PSRSYLKRNPKDFHGYEPEYFLKHQRAHSIEKNCVCSECGKAFRCKSQLIVHLRIHTGERPYECSKCERA
FSAKSNLNAHQVRVHTGEKPYSCSECEKVSFRSQLIVHQEIHGTGKPYGCSECGKAYSWSQQLLHQRS
TGVPKPYECSECGKAFSLKSPFVVHQRTHTGVPKPHKCECGKAFRSKSYLLVHIRMHTGEKPYQCSDCGKA
FNMKTQLIVHQGVHTGNNPYQCGECGKAFRKEQLTAHLRAHAGEKPYGCSECGKAFSSKSYLVIHRRTH
TGERPYECSLCERAFCGKSQLIIHQRTHTSTEKPYECNECEKAYPRKASLQIHQKTHSGEKPFCSECGKA
FTQKSSLSEHQVRVHTGEKPKWCSECGKSFCWNSGLRIHRKTHK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_062537</u>
RefSeq Size:	3198
RefSeq ORF:	1599



[View online »](#)

Synonyms:	HEL-179; KOX20
Locus ID:	7574
UniProt ID:	P17031 , V9HW07
Cytogenetics:	12q24.33
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Protein Families:	Transcription Factors

Product images:



Coomassie blue staining of purified ZNF26 protein (Cat# [TP308543]). The protein was produced from HEK293T cells transfected with ZNF26 cDNA clone (Cat# [RC208543]) using MegaTran 2.0 (Cat# [TT210002]).