

Product datasheet for PH308835


C5orf24 (NM_152409) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	C5orf24 MS Standard C13 and N15-labeled recombinant protein (NP_689622)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC208835
Predicted MW:	20.1 kDa
Protein Sequence:	>RC208835 protein sequence Red=Cloning site Green=Tags(s) MMHPVASSNPAFCGPGKPSCLNEDAMRAADQFDIYSSQQSKYSHTVNHKPMVCQRQDPLNETHLQTTSGR SIEIKDELKKKKLNLRSGKRGRPSGTTKSAGYRTSTGRPLGTTKAAGFKTSPGRPLGTTKAAGYKVPSPGR PPGSIKALSRLADLGYGCGTAAFPYPMMHGRAVHGVEETSSEVKPPNE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_689622
RefSeq Size:	4896
RefSeq ORF:	564
Locus ID:	134553
UniProt ID:	Q7Z6I8
Cytogenetics:	5q31.1



[View online »](#)

Product images:

Coomassie blue staining of purified C5orf24 protein (Cat# [TP308835]). The protein was produced from HEK293T cells transfected with C5orf24 cDNA clone (Cat# [RC208835]) using MegaTran 2.0 (Cat# [TT210002]).