

Product datasheet for PH309529

WWP2 (NM_199424) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	WWP2 MS Standard C13 and N15-labeled recombinant protein (NP_955456)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC209529
Predicted MW:	98.9 kDa
Protein Sequence:	>RC209529 protein sequence Red=Cloning site Green=Tags(s)

MASASSRAGVALPFEKSQLTLKVVSAKPKVHNRQPRINSYVEVAVDGLPSETKKTGKRIGSSELLWNEI
IILNVT AQSHLDLKVWSCHTLRNEELGTASVNL SNVLKNNGGKMEMNQLTLNLQTENKGSVVS GGELTIF
LDGPTVDLGNVPNGSALTDGSQLPSRDSSGTAVAPENRHQPSTNCFGGRSRTHRHSASARTTPATGEQ
SPGARSRHRQPVKNSGHSGLANGTVNDEPTTADPEEPSVVGVTSPPAAPLSVTPNPNTTSLPAPATPAE
GEEPSTSGTQQLPAAAQAPDALPAGWEQREL PNGRVVYVDHNTKTTTWERPLPPGWEKRTDPRGRFYVD
HNTRTTTWRPTAEYVRNYEQWQSQRNQLQGAMQHF SQRFLYQSSASTDHDPLGPLPPGWEKRQDNGRV
YYVNHNRTRTTQWEDPRTQGMIEPALPPGWEMKYTSEGVRYFVDHNRTRTTFKDPRPGFESGKQGS
YDRSFRWKYHQFRFLCHSNALPSHVKISVSRQTLFEDSFQQIMNMKPYDLRRRLYIIMRGEGLDYGGIA
REWFLLSHEVLNPMYCLFEYAGKNYCLQINPASSINPDHLTYFRF IGRFIAMALYHGKFIDTGF TLPF
YKRLNKRPTLKDLESIDPEFYNSIVWIKENNL EECGLELYFIQDMEILGKVTTHELKEGGESIRVTEEN
KEEYIMLLTDWRFTRGV EEQTKAFLDGFNEVAPLEWLR YFDEKELEMLCGMQEIDMSDWQKSTIYRHYT
KNSKQIQWFQVVKEMDNEKRIRLLQFVTGTCRLPVGGFAELIGSNGPQKFCIDKVGKETWLP RSHTCFN
RLDLPPYKSYEQLREKLLYAI EETEGFGQE

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.



[View online »](#)

RefSeq: [NP_955456](#)

RefSeq Size: 3700

RefSeq ORF: 2613

Synonyms: AIP2; WWp2-like

Locus ID: 11060

UniProt ID: [O00308](#), [B4DHF6](#)


Cytogenetics: 16q22.1

Summary: This gene encodes a member of the Nedd4 family of E3 ligases, which play an important role in protein ubiquitination. The encoded protein contains four WW domains and may play a role in multiple processes including chondrogenesis and the regulation of oncogenic signaling pathways via interactions with Smad proteins and the tumor suppressor PTEN. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 10. [provided by RefSeq, Jul 2012]

Protein Families: Druggable Genome

Protein Pathways: Ubiquitin mediated proteolysis

Product images:



Coomassie blue staining of purified WWP2 protein (Cat# [TP309529]). The protein was produced from HEK293T cells transfected with WWP2 cDNA clone (Cat# [RC209529]) using MegaTran 2.0 (Cat# [TT210002]).