

## Product datasheet for PH309704

### NECAB1 (NM\_022351) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	NECAB1 MS Standard C13 and N15-labeled recombinant protein (NP_071746)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC209704
Predicted MW:	40.5 kDa
Protein Sequence:	>RC209704 protein sequence Red=Cloning site Green=Tags(s)
	MEDSQETSPSSNNSSEELSSALHLKSGMSIFLDILRRADKNDGKLSFEEFKAYFADGVLSGEELHELFB TIDHTNNTLDEELCEYFSQHLGEYENVLALEGLNLKAMGKTKKDYQEASNLEQFVTRFLLKETL NQLQSLQNSLECAMETTEEQTRQERQGPAPKPEVLSIQWPGKRSSRRVQRHNSFSPNSPQFNVSQGLLEE DNQWMTQINRLQKLDRLKDLKLEPPEEEIEGNTKSHIMLVQRQMSVIEEDLEEFQLALKHYVESAS SQSGCLRISIQKLSNESRYMIYEFWENSSVWNSHLQTNYSKTFQRSNVDFLETPELTSTMLVPASWWILN N
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>4</sub> ]-L-Arginine and [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>2</sub> ]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u><a href="#">NP_071746</a></u>
RefSeq Size:	5140
RefSeq ORF:	1053
Synonyms:	EFCBP1; STIP-1
Locus ID:	64168




[View online »](#)

UniProt ID: [Q8N987](#)

Cytogenetics: 8q21.3

### Product images:



Coomassie blue staining of purified NECAB1 protein (Cat# [TP309704]). The protein was produced from HEK293T cells transfected with NECAB1 cDNA clone (Cat# [RC209704]) using MegaTran 2.0 (Cat# [TT210002]).