

Product datasheet for PH309795

67kDa Laminin Receptor (RPSA) (NM_002295) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	RPSA MS Standard C13 and N15-labeled recombinant protein (NP_002286)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC209795
Predicted MW:	32.8 kDa
Protein Sequence:	>RC209795 protein sequence Red=Cloning site Green=Tags(s) MSGALDVLQMKEEDVLKFLAAGTHLGGTNLDFQMEQYIYKRKSDGIYIINLKRTWEKLLLAAARIVAIEN PADVSVISSRNTGQRAVLKFAAATGATPIAGRFTPGFTTNQIQAAFREPRLLVVTDPRADHQPLTEASYV NLPTIALCNTDSPLRYVDIAIPCNNKGAHSVGLMWWMLAREVLRMRGTISREHPWEVMPDLYFYRDPPEEI EKEGQAAA EKAVTKKEEFQGEWTAPAPEFTATQPEVADWSEGVQVPSVPIQQFPTEDWSAQPATEDWSAAP TAQATEWVGATTDWS TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_002286
RefSeq Size:	1155
RefSeq ORF:	885
Synonyms:	37LRP; 67LR; ICAS; LAMBR; lamR; LAMR1; LBP; LBP/p40; LRP; LRP/LR; NEM/1CHD4; p40; SA
Locus ID:	3921



[View online »](#)

UniProt ID: [P08865](#), [A0A024R2L6](#)


Cytogenetics: 3p22.1

Summary: Laminins, a family of extracellular matrix glycoproteins, are the major noncollagenous constituent of basement membranes. They have been implicated in a wide variety of biological processes including cell adhesion, differentiation, migration, signaling, neurite outgrowth and metastasis. Many of the effects of laminin are mediated through interactions with cell surface receptors. These receptors include members of the integrin family, as well as non-integrin laminin-binding proteins. This gene encodes a high-affinity, non-integrin family, laminin receptor 1. This receptor has been variously called 67 kD laminin receptor, 37 kD laminin receptor precursor (37LRP) and p40 ribosome-associated protein. The amino acid sequence of laminin receptor 1 is highly conserved through evolution, suggesting a key biological function. It has been observed that the level of the laminin receptor transcript is higher in colon carcinoma tissue and lung cancer cell line than their normal counterparts. Also, there is a correlation between the upregulation of this polypeptide in cancer cells and their invasive and metastatic phenotype. Multiple copies of this gene exist, however, most of them are pseudogenes thought to have arisen from retropositional events. Two alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Ribosome

Product images:



Coomassie blue staining of purified RPSA protein (Cat# [TP309795]). The protein was produced from HEK293T cells transfected with RPSA cDNA clone (Cat# [RC209795]) using MegaTran 2.0 (Cat# [TT210002]).