


Product datasheet for PH311067

TXLNGY (NM_001005852) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	CYorf15A MS Standard C13 and N15-labeled recombinant protein (NP_001005852)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC211067
Predicted MW:	14.7 kDa
Protein Sequence:	>RC211067 protein sequence Red=Cloning site Green=Tags(s) MEEAGLCGLREKADMLCNSESHDILQHQDSNCSATSNKHLLDEEGRDFITKNRSWVSPVHCTQESRREL PEQEVAPPSGQQALQCNRNKEKVLERQGLSLSPRLERGGTFMGHGSCLKPELVILSLQTLE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001005852
RefSeq Size:	872
RefSeq ORF:	393
Synonyms:	FLJ33216; MGC131732
Locus ID:	246126
Cytogenetics:	Yq11.222-q11.223


[View online »](#)

Product images:


Coomassie blue staining of purified CYorf15A protein (Cat# [TP311067]). The protein was produced from HEK293T cells transfected with CYorf15A cDNA clone (Cat# [RC211067]) using MegaTran 2.0 (Cat# [TT210002]).