

Product datasheet for PH311269

C20orf79 (SCP2D1) (NM_178483) Human Mass Spec Standard

Product data:


Product Type:	Mass Spec Standards
Description:	C20orf79 MS Standard C13 and N15-labeled recombinant protein (NP_848578)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC211269
Predicted MW:	17.7 kDa
Protein Sequence:	>RC211269 protein sequence Red=Cloning site Green=Tags(s) MWKRSDHQPKIKAEDGPLVGQFEVLGSPPEPAMPHLELSEFESFPVFQDIRLHIREVGAQLVKKVNAVF QLDITKNGKTIILRWTIDLKNGSGMYPGPARLPADTVFTIPESVFMELVLGKMNPQKAFLAGKFKVSGKV LLSWKLERVFKDWAKF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_848578
RefSeq Size:	677
RefSeq ORF:	468
Synonyms:	C20orf79; dj1068E13.2; HSD22
Locus ID:	140856
UniProt ID:	Q9UJQ7



[View online »](#)

Cytogenetics: 20p11.23

Product images:



Coomassie blue staining of purified SCP2D1 protein (Cat# [TP311269]). The protein was produced from HEK293T cells transfected with SCP2D1 cDNA clone (Cat# [RC211269]) using MegaTran 2.0 (Cat# [TT210002]).