

Product datasheet for PH311599

CIAP2 (BIRC3) (NM_182962) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	BIRC3 MS Standard C13 and N15-labeled recombinant protein (NP_892007)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC211599
Predicted MW:	68.4 kDa
Protein Sequence:	>RC211599 protein sequence Red=Cloning site Green=Tags(s)

MNIVENSIFLSNLMKSANTFELKYDLSCELYRMSTYSTFPAGVPVRSERSLARAGFYTGVDKVKCFCCG
LMLDNWKRGDSPTEKHKLYPSCRFFVQSLNSVNNLEATSQPTFPSSVTNSTHLLPGTENSIFYFRGSYSN
SPSNPVNSRANQDF SALMRSSYHCAMNNENARLLTFQTWPLTFLSPTDLAKAGFYIIGPGDRVACFACGG
KLSNWEPKDNAMSEHLRHFPKCPF IENQLQDTSRYTVSNLSMQTHAARFKTFNWPSSVLVNPEQLASAG
FYYVGNSSDDVKCFCCDGGRLRCWESGDDPWVQHAKWFPRCEYLIRIKGQEFIRQVQASYPHLLEQLLSTSD
SPGDENAESSIIHFEPGEDHSEDAIMMNTPVINA AAVEMGFSRSLVKQTVQRKILATGENYRLVNDLVLDL
LNAEDEIREEERERATEEKESENDLLLIRKNRMALFQHLTCVIPILDSLLTAGIINEQEHVDIKQKTQTSL
QARELIDTILVKGNI AATVFRNSLQEA EAVLYEHLFVQQDIKIYPTEDVSDLPVEEQLRRLQEERTCKVC
MDKEVSIVFIPCGHLVCKDCAPSLRKCPICRSTIKGTVRTFLS


TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_892007
RefSeq Size:	4372
RefSeq ORF:	1812



[View online »](#)

Synonyms:	AIP1; API2; c-IAP2; CIAP2; HAIP1; HIAP1; IAP-1; MALT2; MIHC; RNF49
Locus ID:	330
UniProt ID:	Q13489
Cytogenetics:	11q22.2
Summary:	This gene encodes a member of the IAP family of proteins that inhibit apoptosis by binding to tumor necrosis factor receptor-associated factors TRAF1 and TRAF2, probably by interfering with activation of ICE-like proteases. The encoded protein inhibits apoptosis induced by serum deprivation but does not affect apoptosis resulting from exposure to menadione, a potent inducer of free radicals. It contains 3 baculovirus IAP repeats and a ring finger domain. Transcript variants encoding the same isoform have been identified. [provided by RefSeq, Aug 2011]
Protein Families:	Druggable Genome
Protein Pathways:	Apoptosis, Focal adhesion, NOD-like receptor signaling pathway, Pathways in cancer, Small cell lung cancer, Ubiquitin mediated proteolysis

Product images:

Coomassie blue staining of purified BIRC3 protein (Cat# [TP311599]). The protein was produced from HEK293T cells transfected with BIRC3 cDNA clone (Cat# [RC211599]) using MegaTran 2.0 (Cat# [TT210002]).