

Product datasheet for PH312273

JNK2 (MAPK9) (NM_139070) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	MAPK9 MS Standard C13 and N15-labeled recombinant protein (NP_620709)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC212273
Predicted MW:	48.1 kDa
Protein Sequence:	>RC212273 representing NM_139070 Red=Cloning site Green=Tags(s)
	MSDSKCDSQFYYSVQVADSTFTVLKRYQQLKPIGSGAQGIVCAAFDVLGINVAVKLSRPFQNTAKRA YRELVLLKCVNHKNIISLLNVFTPQKTLLEEFQDVYLVMEMLDANLCQVIHMELDHERMSYLLYQMLCGIK HLHSAGIIHRDLKPSNIVVKSCTLKILDFGLARTACTNFMTPYVVTRYRAPEVILGMGYKENVDIWS VGCIMAEMVLHKVLPGRDYIDQWNVIEQLGTPSAEFMKKLQPTVRNYVENRPKYPGIKFEELFPDWIF PSESERDKIKTSQARDLLSKMLVIDPKRISVDEALRHPYITVWYDPAEAEAPPPQIYDAQLEEREHAIE EWKELIYKEVMDWEERSKNGVVKDQPSDAAVSSNATPSQSSSINDISSMSTEQTLASDTDSSLDASTGPL EGCR
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_620709</u>
RefSeq Size:	1942
RefSeq ORF:	1272
Synonyms:	JNK-55; JNK2; JNK2A; JNK2ALPHA; JNK2B; JNK2BETA; p54a; p54aSAPK; PRKM9; SAPK; SAPK1a



[View online »](#)

Locus ID: 5601

UniProt ID: [P45984](#)


Cytogenetics: 5q35.3

Summary: The protein encoded by this gene is a member of the MAP kinase family. MAP kinases act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. This kinase targets specific transcription factors, and thus mediates immediate-early gene expression in response to various cell stimuli. It is most closely related to MAPK8, both of which are involved in UV radiation induced apoptosis, thought to be related to the cytochrome c-mediated cell death pathway. This gene and MAPK8 are also known as c-Jun N-terminal kinases. This kinase blocks the ubiquitination of tumor suppressor p53, and thus it increases the stability of p53 in nonstressed cells. Studies of this gene's mouse counterpart suggest a key role in T-cell differentiation. Several alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Sep 2008]

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase

Protein Pathways: Adipocytokine signaling pathway, Colorectal cancer, Epithelial cell signaling in Helicobacter pylori infection, ErbB signaling pathway, Fc epsilon RI signaling pathway, Focal adhesion, GnRH signaling pathway, Insulin signaling pathway, MAPK signaling pathway, Neurotrophin signaling pathway, NOD-like receptor signaling pathway, Pancreatic cancer, Pathways in cancer, Progesterone-mediated oocyte maturation, RIG-I-like receptor signaling pathway, T cell receptor signaling pathway, Toll-like receptor signaling pathway, Type II diabetes mellitus, Wnt signaling pathway

Product images:



Coomassie blue staining of purified MAPK9 protein (Cat# [TP312273]). The protein was produced from HEK293T cells transfected with MAPK9 cDNA clone (Cat# [RC212273]) using MegaTran 2.0 (Cat# [TT210002]).