

Product datasheet for PH312450

C20orf151 (RBBP8NL) (NM_080833) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	C20orf151 MS Standard C13 and N15-labeled recombinant protein (NP_543023)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC212450
Predicted MW:	71.3 kDa
Protein Sequence:	>RC212450 representing NM_080833 Red=Cloning site Green=Tags(s)

MESFMESLNRLKEIHEKEVLGLQNKLELNSERCDAQRIEELFSKNHQLREQQKTLKENLRVLENRLRA
GLCDRCMVTQELARKRQQEFESSHLQNLQRIFILTNEMNGLKEENETLKEEVKRLRGLGDRPKPRAKEGT
SDPPSPLLLPSPPGGWKAITEKPPGGHEEAEEDHQVGLRGEEKPAGHRTSPVAKISPGATLPESRAPDMS
PQRISNQLHGTIAVVRPGSQACPADRGANGTPPPLPARSSPPSPAYERGLSLDSFLRASRPSAMTHEAP
KLSPKVDRLCLLNRPLSLHLQSPHSSPLAPAAAPSDPRLQDLKAREAEAWEEPTELLGLPSALAGMQDLR
LEGALHLLLAQQQLRARARAGSVRPRGQPTPGEMLPSPVGSDEGPENEGTRAALAAAGLSGGRHTQPA
GPGRAQRTEAAATQDCALDKPLDLSEWGRARGQDTPKAGQHGSLSAAAHTASPEPPTQSGPLTRSPQA
LSNGTKGTRVPEQEEASTPMDPSRPLPGSQLSLSSPGSTEDEDTGRPLPPHPQPPHPQPPDLGHPEP
SKAEVLRPESELDDETDTPGSEVLSSQAEATTSTTGEGPECICTQEHGQGPPrKRKRASEPGDKASKKP
SRGRRKL TATEGPSRDAEDHSPSPNSSPWEET

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_543023
RefSeq Size:	2793



[View online »](#)

RefSeq ORF:	1992
Synonyms:	C20orf151
Locus ID:	140893
UniProt ID:	Q8NC74
Cytogenetics:	20q13.33

Product images:



Coomassie blue staining of purified RBBP8NL protein (Cat# [TP312450]). The protein was produced from HEK293T cells transfected with RBBP8NL cDNA clone (Cat# [RC212450]) using MegaTran 2.0 (Cat# [TT210002]).