

Product datasheet for PH312603


FAM177B (NM_207468) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	FAM177B MS Standard C13 and N15-labeled recombinant protein (NP_997351)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC212603
Predicted MW:	18 kDa
Protein Sequence:	>RC212603 representing NM_207468 Red=Cloning site Green=Tags(s) MEIDGFQQLDLEKSVPSKTTTPKRIIHFDVGDIMEEYSTEETEEEEKEEQSTNSTLDPSKLSWGPYLRFW AGRIASTSFSTCEFLGGRFAVFFGLTQPKYQYVLNEFYRIQNKSDNKSERRGSKAQAAEVPNEKCHLEA GVQEYGTIQQDVTEAIPQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_997351</u>
RefSeq Size:	1345
RefSeq ORF:	474
Locus ID:	400823
UniProt ID:	<u>A6PVY3</u>
Cytogenetics:	1q41


[View online »](#)

Product images:



Coomassie blue staining of purified FAM177B protein (Cat# [TP312603]). The protein was produced from HEK293T cells transfected with FAM177B cDNA clone (Cat# [RC212603]) using MegaTran 2.0 (Cat# [TT210002]).