

Product datasheet for PH313161

ARSA (NM_001085427) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	ARSA MS Standard C13 and N15-labeled recombinant protein (NP_001078896)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC213161
Predicted MW:	53.6 kDa
Protein Sequence:	>RC213161 protein sequence Red=Cloning site Green=Tags(s)

MGAPRSLLLALAAGLAVARPPNIVLIFADDLGYGDLGCGYHPSSTTPNLDQLAAGGLRFTDFYVPVSLCT
PSRAALLTGRLPVRMGMPGVLPSSRGGPLEEVTVAEVLAAARGYL TGMAGKWHLGVGPEGAFLPPHQG
FHRFLGIPYSHDQGPCQNLTCFPPATPCDGGCDQGLVPIPLLANLSVEAQPWLPGLEARYMAFAHDLMA
DAQRQDRPFLLYASHHHTYPQFSGQSF AERSGRGPFGDSLMELDAAVGTLMTAIGDLGLEETLVIFTA
DNGPETMRMSRGGCSGLLRGCGTTEYEGVREPALAFWPGHIAPGVTHELASSLDLLPTLAALAGAPLPN
VTLDGFDLSPLLLGTGKSPRQSLFFYPSYPDEVRGVFAVRSGKYKAHFF TQGSAHSDTTADPACHASSSL
TAHEPPLLYDL SKDPGENYNLLGGVAGATPEVLQALKQLQLLKAQLDAAVTFGPSQVARGEDPALQICCH
PGCTPRPACCHCPDPA

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001078896
RefSeq Size:	4127
RefSeq ORF:	1521



[View online »](#)

Synonyms:	ASA; MLD
Locus ID:	410
UniProt ID:	P15289 , A0A0C4DFZ2
Cytogenetics:	22q13.33
Summary:	The protein encoded by this gene hydrolyzes cerebroside sulfate to cerebroside and sulfate. Defects in this gene lead to metachromatic leucodystrophy (MLD), a progressive demyelination disease which results in a variety of neurological symptoms and ultimately death. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Dec 2010]
Protein Families:	Druggable Genome
Protein Pathways:	Lysosome, Sphingolipid metabolism

Product images:



Coomassie blue staining of purified ARSA protein (Cat# [TP313161]). The protein was produced from HEK293T cells transfected with ARSA cDNA clone (Cat# [RC213161]) using MegaTran 2.0 (Cat# [TT210002]).