

## Product datasheet for PH313633

### TMEM98 (NM\_001033504) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	TMEM98 MS Standard C13 and N15-labeled recombinant protein (NP_001028676)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC213633
Predicted MW:	24.6 kDa
Protein Sequence:	>RC213633 protein sequence Red=Cloning site Green=Tags(s)  METVVIVAIGVLATIFLASFAALVLVCRQRYCRPRDLLQRYDSKPIVDLIGAMETQSEPSELELDDVVIT NPHIEAILENEDWIEDASGLMSHCIAILKICHTLTEKLVAMTMGSGAKMKT SASVSDIIVVAKRISPRVD DVYKSMYPPLDPKLLDARTTALLLSVSHLVLVTRNACHLTGGLDWIDQSL SAAEEHLEVLREAALASEPD KGLPGPEGFLQEQSAI  TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>4</sub> ]-L-Arginine and [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>2</sub> ]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<a href="#">NP_001028676</a>
RefSeq Size:	1732
RefSeq ORF:	678
Synonyms:	TADA1
Locus ID:	26022
UniProt ID:	<a href="#">Q9Y2Y6</a>




[View online »](#)

Cytogenetics: 17q11.2

Summary: This gene encodes a transmembrane protein. A missense mutation in this gene result in Nanophthalmos 4 (NNO4). Alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Sep 2014]

Protein Families: Transmembrane

### Product images:



Coomassie blue staining of purified TMEM98 protein (Cat# [TP313633]). The protein was produced from HEK293T cells transfected with TMEM98 cDNA clone (Cat# [RC213633]) using MegaTran 2.0 (Cat# [TT210002]).