

## Product datasheet for PH313966

### NFAT1 (NFATC2) (NM\_173091) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	NFATC2 MS Standard C13 and N15-labeled recombinant protein (NP_775114)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC213966
Predicted MW:	100 kDa
Protein Sequence:	>RC213966 representing NM_173091 Red=Cloning site Green=Tags(s)

MNAPERQPQPDGGDAPGHEPGGSPQDELDFSIILFDYEYLNPNEEEEPNAHKVASPPSPGAYPDDVLDYGLK  
PYSPLASLSGEPGRFGEPDRVGPQKFLSAAKPAGASGLSPRIEITPSHELIQAVGPLRMRDAGLLVEQP  
PLAGVAASPRFTLPVPGFEGYREPLCLSPASSGSSASFISDTFSPYTPSPCVSPNNGPDDLCQFQNI  
HYSPTSPIMSPRTSLAEDSCLGRHSPVPRPASRSSSPGAKRRHSCEALVALPPGASPQRSRSPSPQPS  
SHVAPQDHGSPAGYPPVAGSAVIMDALNSLATDSPCGIPPKMWTSPDPSVSAAPSKAGLPRHIYPAVE  
FLGPCEQGERRNSAPESILLVPPTWPKPLVPAIPICSIPVTASLPLEWPLSSQSGSYELRIEVQPKPHH  
RAHYETEGSRGAVKAPTGGHPVVQLHGYMENKPLGLQIFIGTADERILKPHAFYQVHRITGKTVTTTSYE  
KIVGNTKVL EIPLEPKNNMRATIDCAGILKLRNADIELRKGETDIGRKNTRVRLVFRVHIPESSGRIVSL  
QTASNPIECSQRSHELPMVERQDTSCLVYGGQMILTGQNF TSESKVVFTEKTTDGGQIWEMEATVDK  
DKSQPNMLFVEIPEYRNKHIRTVPKVN FYINGKRKRSQPQHFTYHPVPAIKTEPTDEYDPTLICSPTHG  
GLGSQPYYPQHMPVAESPSCLVATMAPCQFRTGLSSPDARYQQNPAAVLYQRSKSLSPSLLGYQQPAL  
MAAPLSLADAHRSVLVHAGSQGQSSALLHPSPTNQQASPVIIHYSPTNQQLRGSHQEFQHIMYCENFAPG  
TTRPGPPPVSQGQRLSPGSYPTVIQQNATSQRAAKNGPPVSDQKEVLPAGVTIKQEQLNDQTYLDDVNE  
IIRKEFSGPPARNQT

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>4</sub> ]-L-Arginine and [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>2</sub> ]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.



[View online »](#)

RefSeq: [NP\\_775114](#)

RefSeq Size: 3166

RefSeq ORF: 2775

Synonyms: NFAT1; NFATP

Locus ID: 4773

UniProt ID: [Q13469](#)


Cytogenetics: 20q13.2

**Summary:** This gene is a member of the nuclear factor of activated T cells (NFAT) family. The product of this gene is a DNA-binding protein with a REL-homology region (RHR) and an NFAT-homology region (NHR). This protein is present in the cytosol and only translocates to the nucleus upon T cell receptor (TCR) stimulation, where it becomes a member of the nuclear factors of activated T cells transcription complex. This complex plays a central role in inducing gene transcription during the immune response. Alternate transcriptional splice variants encoding different isoforms have been characterized. [provided by RefSeq, Apr 2012]

**Protein Families:** Druggable Genome, Transcription Factors

**Protein Pathways:** Axon guidance, B cell receptor signaling pathway, MAPK signaling pathway, Natural killer cell mediated cytotoxicity, T cell receptor signaling pathway, VEGF signaling pathway, Wnt signaling pathway

### Product images:



Coomassie blue staining of purified NFATC2 protein (Cat# [TP313966]). The protein was produced from HEK293T cells transfected with NFATC2 cDNA clone (Cat# [RC213966]) using MegaTran 2.0 (Cat# [TT210002]).